


RAPid® Zap


Warnings and Instructions: Review carefully, proper application is required

⚠ Warning: This device will not prevent or eliminate risk of injury. Do Not Overtighten. If swelling, pain, skin irritation, or an unusual reaction occurs, discontinue use immediately and consult your medical professional.


Care: Hand wash using mild soap. Rinse thoroughly. Air dry only. Do not tumble dry.




1. This brace is designed to be worn over a sock. Loosen all straps and laces and place foot in brace.




2. Position your heel snugly at the rear bottom corner of the brace. Pull and fasten two lower straps.




3. Fasten one end of the upper strap. Pull and fasten other end of upper strap.




4. Adjust straps as necessary for a snug fit.




5. Bring the inside strap forward, across top of foot, then under the heel.




6. Pull and secure strap on inside of ankle.




7. Bring the outside strap forward, across top of foot, then under the heel. Pull and secure strap on outside of ankle.



8. Stand, putting full weight on ankle. Pull up on both straps simultaneously and secure.



9. Pull upper elastic strap from rear and fasten in front.



10. Completed application.



Scan for video

HELY & WEBER
2023 Hely & Weber

