



Warnings and Instructions: Review carefully, proper application is required

⚠ Warning: This device will not prevent or eliminate risk of injury. Do Not Overtighten. If swelling, pain, skin irritation, or an unusual reaction occurs, discontinue use immediately and consult your medical professional.


Care: Hand wash using mild soap. Rinse thoroughly. Air dry only. Do not tumble dry. For single patient use only.




1. Center brace over patella.




2. Pull lower strap around back of leg, through loop lock, and pull and fasten.




3. Pull upper straps around back of leg, through loop locks, and pull and fasten.




4. Undo both outer flaps. Make sure hinges are aligned with knee joint.




5. Place the buttress in the desired position.




6. Hold buttress in place & pull outer flap straps to desired tension.



7. Fasten outer flap straps.



8. Pull and fasten inner straps.



9. Finished application.

