


Shoulder Houlder II Universal


Warnings and Instructions: Review carefully, proper application is required

⚠ Warning: This device will not prevent or eliminate risk of injury. Do Not Overtighten. If swelling, pain, skin irritation, or an unusual reaction occurs, discontinue use immediately and consult your medical professional.


Care: Hand wash using mild soap. Rinse thoroughly. Air dry only. Do not tumble dry.




1. Adjust sling length to length of olecranon to knuckles by folding and fastening the front edge of the sling.




2. Place pillow against body and secure waist strap around waist.




3. Place affected elbow in corner of elbow pocket, and attach pocket to pillow.




4. Fasten ends of shoulder strap to front corners of brace.
Secure arm in brace with elbow strap



5. Have practitioner adjust back straps before first time use.



6. Make final adjustments to brace using front straps.
Secure thumb strap.



7. Completed application. Remove brace using quick release buckles.



Scan to view application video

