

Webly® Achilles TT Application

Warnings and Instructions: Review carefully, proper application is required


Warning: This device will not prevent or eliminate risk of injury. Do Not Overtighten. If swelling, pain, skin irritation, or an unusual reaction occurs, discontinue use immediately and consult your medical professional. For optimal performance wear over a sock.
Care: Hand wash using mild soap. Rinse thoroughly. Air dry only. Do not tumble dry.




1. Preparation for Fitting: Pull BOA® dial to unlock. Pull on both sides of gauntlet to loosen straps. Remove posterior panel.




2. Slide foot into gauntlet with laces loosely tightened.




3. Center the Achilles tendon between the two soft pads and attach one side of the posterior panel to the ankle gauntlet and then the other.




4. When back panel is secured, push in and turn the BOA® dial clockwise to tighten the lacing to create intimate fit of the gauntlet. Note: When tightened, pressure on both sides of the Achilles tendon should be felt.




5. Pull the tendon tension adjustment tab to a comfortable tension so it compresses the tendon between the two pads (do not overtighten).



6. Application is complete.



7. To remove the Webly® Achilles TT: Loosen tendon strap, pull out on the BOA® dial, open the laces entirely.



8. Undo lateral or medial side of the back panel and slide foot out OR slide brace over foot.

HELY&WEBER

2025 Hely & Weber

To Reapply: Slide foot in brace, make sure the two posterior pads approximate the tendon, tighten the BOA®, cinch the tendon tension adjustment tab to a comfortable tension.

Application Video