GK ORIGIN

A SENSIBLE SOLUTION FOR CONTROLLING PAIN OF OA PATIENTS



5640 GK ORIGIN

INDICATIONS:

- MCL/LCL INSTABILITY
- PATELLOFEMORAL TRACKING MALALIGNMENT
- MULTI-COMPARTMENT OSTEOARTHRITIS
- CHONDROMALACIA PATELLA

FEATURES

- PATELLA SLING LIFTS AND PROVIDES DIRECTIONAL CONTROL OF THE PATELLA REDUCING DISCOMFORT
- ADJUSTABLE CALF ALLOWS THE BRACE TO ACCOMMODATE EVEN THE HARD TO FIT PATIENTS
- KUHL NEOPRENE WICKS AWAY MOISTURE WHILE PROVIDING WARMTH AND COMPRESSION
- BUILT IN CONDYLE PADS IMPROVE COMFORT AND INCREASE PATIENT COMPLIANCE
- MEETS MEDICARE'S CONSERVATIVE CARE PROTOCOL PRIOR TO TKR

A SOLUTION WHEN INSTABILITY CAN'T BE DOCUMENTED





🔘 in



PDAC APPROVER

GK ORIGIN

PROVIDES SUPPORT, STABILITY, COMPRESSION AND PAIN RELIEF

5640 GK ORIGIN

The patella sling and "U" buttress lift and control the patella with an inferior to superior directional pull. This places the patella in a more natural position and helps patients move with more stability and less pain.

Fits Real Patients!

The adjustable calf allows the brace to be fit on a conical shaped leg. Accommodating a wide range of patients is now easier.





BREATHABLE NEOPRENE AND BUILT IN CONDYLE PAD

ORDERING INFORMATION

ORDER#	SIZE	MIN IN	MAX IN	MIN CM	MAX CM
5640	XS	13″	15″	33CM	38CM
5640	S	15″	17″	38CM	43CM
5640	Μ	17″	20″	43CM	51CM
5640	L	19″	22″	48CM	56CM
5640	XL	21″	24″	53CM	61CM
5640	XXL	23″	27″	58CM	69CM
DDAC A.		14000			

PDAC Approved L1820 THIGH: MEASURE 6 INCHES ABOVE PATELLA

CALL TODAY! (800) 654-3241 EMAIL: ORDER@HELY-WEBER.COM



© 2018 HELY & WEBE