

NoSwet[®] Knapp[™]

HINGED KNEE ORTHOSIS



A NEW SOLUTION FOR KNEE STABILIZATION

7676/7678 NOSWET[®] KNAPP[™]

INDICATIONS:

- STRAINS AND SPRAINS OF THE MCL/LCL, ACL/PCL
- MENISCUS INJURIES

FEATURES AND BENEFITS:

- **HYBRID OF DOUBLE KNIT AND NOSWET[®] MATERIALS** - ALLOWS FOR GREAT BREATHABILITY AND COMFORT.
- **LYCRA BACK PANEL** RESTORES COMPRESSION TO NON-NEOPRENE BRACES TO REDUCE SWELLING AND KEEP THE ATHLETE COOL.
- **ANTERIOR CLOSURE DESIGN** - PROVIDES EASY ON AND OFF AND A CUSTOMIZED FIT.
- **INTERNAL HINGES** - FOR ATHLETIC EVENTS THAT REQUIRE COVERED HINGES, MOLDABLE ALUMINUM HINGES ARE INSIDE POCKETS SO THERE IS NO EXPOSED METAL.
- **ELASTICIZED STRAPS** - ALLOWS THE STRAPS TO MOVE WITH THE MUSCLES PROVIDING AN INTIMATE FIT AS WELL AS ENHANCED SUSPENSION.
- **TWO AVAILABLE LENGTHS** - 12" AND 16" LENGTHS PROVIDE OPTIMAL LIGAMENOUS SUPPORT.

DESIGNED IN
COLLABORATION WITH:



THOMAS P. KNAPP, MD
CEDARS-SINAI KERLAN-JOBE INSTITUTE

HELY & WEBER
BRACING INNOVATIONS



**DOUBLE KNIT
BACK PANEL
AND BUILT IN
CONDYLE PADS
PROVIDE
COMFORT,
SUPPORT AND
UNINHIBITED
MOBILITY.**



**SLEEK DESIGN PROVIDES A
COMFORTABLE FIT WITH
COVERED HINGES AND A
SPACERMESH BODY.**

**ANTERIOR
CLOSURE
DESIGN
PROVIDES A
CUSTOMIZED
FIT.**



**TWO LENGTHS
AND TWO HINGE
OPTIONS.
12" AND 16"
LENGTHS,
STANDARD OR
HEAVY-DUTY
HINGES FOR
COMPLETE
VERSATILITY
AND USE.**



7676/7678 ORDERING INFORMATION:

SIZING CHART		
XS	13"-15"	(33CM - 38CM)
S	15"-17"	(38CM - 43CM)
M	17"-19"	(43CM - 48CM)
L	19"-22"	(48CM - 56CM)
XL	22"-25"	(56CM - 63CM)
XXL	25"-28"	(63CM - 71CM)

PART #7678 IS FOR THE LONG VERSION OF THE BRACE.
MEASURE AT 6" ABOVE MID-PATELLA WITH LEG FULLY EXTENDED.

CALL TODAY! (800) 654-3241
EMAIL: ORDER@HEL-WEBER.COM

