

# Pat-Trac™

## HINGED KNEE ORTHOSIS



GET **THE** SUPPORT YOUR KNEE  
DESERVES

### 5672 PAT-TRAC™

#### INDICATIONS:

- PATELLA SUBLUXATION
- ANTERIOR KNEE PAIN
- PATELLAR MALTRACKING
- CHONDROMALACIA PATELLA
- PATELLA BAJA

#### FEATURES AND BENEFITS:

- **REMOVABLE/ADJUSTABLE PATELLAR STRAP** – ENSURES PATELLAR BUTTRESS IS IN THE OPTIMAL POSITION FOR PATELLAR SUPPORT.
- **THREE UNIVERSAL SIZES** – UNIQUE DESIGN REDUCES SKUS. 3 SIZES AVAILABLE (DESIGNED TO BE USED ON EITHER LEG) REDUCING INVENTORY NEEDS.
- **ANTERIOR CLOSURE DESIGN** – PROVIDES EASY DON AND DOFF AND A CUSTOMIZED FIT.
- **INTERNAL HINGES** – CONFORM TO THE PATIENT'S LEG FOR AN INTIMATE FIT.
- **CONDYLAR PADS** – PROVIDE PADDING OVER HINGE ASSEMBLIES AS WELL AS ENHANCED COMFORT AND SUPPORT AROUND THE KNEE.





LARGE BUILT-IN BUTTRESS MAINTAINS PATELLAR POSITION DURING ACTIVITIES

REMOVABLE PLASTIC HINGE CONFORMS TO THE BODY FOR INTIMATE FIT

ADJUSTABLE PATELLAR FLAP ENSURES OPTIMUM BUTTRESS PLACEMENT

THE PAT-TRAC™ PROVIDES THE EXTRA SUPPORT RECOVERING ATHLETES NEED TO FIND THEIR COMFORT ZONE.



1



2



3



**5672 ORDERING INFORMATION:**

SIZING: MEASUREMENTS ARE TAKEN 4" ABOVE MID-PATELLA		
XS/S	11"-16"	(28CM - 40CM)
M/L	15"-22"	(38CM - 56CM)
XL/XXL	19"-27"	(48CM - 68CM)

CALL TODAY! (800) 654-3241  
EMAIL: [ORDER@HELY-WEBER.COM](mailto:ORDER@HELY-WEBER.COM)

