ORDER NO.

KNAPP™ HINGED KNEE ORTHOSIS - ANTERIOR CLOSURE

KNAPP™ HINGED KNEE ORTHOSIS KUHL™ NEOPREN 5656 5658 12" STANDARD HINGES 16" STANDARD HINGES

12" HEAVY HINGES 16" HEAVY HINGES 5656HH 5658HH

L1820/APPROVED

12" RANGE OF MOTION 5656ROM 16" RANGE OF MOTION 5658ROM L1832/L1833/APPROVED

NOSWET® KNAPP™ HINGED KNEE ORTHOSIS 7656 12" STANDARD HINGES 7658 16" STANDARD HINGES

12" HEAVY HINGES 16" HEAVY HINGES 7656HH 7658HH PDAC L1820 APPROVED

12" RANGE OF MOTION 7656ROM PDAC L1832/L1833 APPROVED

| SIZING: | MEASURE AT MID-PATELLA | WITH LEG FULLY EXTENDED |
|---------|------------------------|-------------------------|
| XS | 10" - 12" | (25CM - 31CM) |
| S | 12" - 14" | (31CM - 36CM) |
| М | 14" - 16" | (36CM - 41CM) |
| L | 16" - 18" | (41CM - 46CM) |
| XL | 18" - 20" | (46CM - 51CM) |
| XXL | 20" - 22" | (51CM - 56CM) |
| | | |

Larger & custom sizes for 5656/5658 MODELS AVAILABLE UPON REQUEST

HOMAS P. KNAPP, M.D.



ORDER NO.

Cedars-Sinai® Kerlan-Jobe Institute



- 3 hinge styles
- Easily accessible hinge allows interchangeability of hinge styles
- Hinges are fully covered for sports
- Elasticized thigh strap ensures comfort
- Ideal for pre/post-operative care of MCL/LCL and meniscal repairs
- Longer length can be used for mild to moderate ACL/PCL sprains & strains and as an adjunct to patella fracture



VELOCITY HINGED KNEE - PULL-UP/ANTERIOR CLOSURE

5656/7656



| | LEG FULLY | EXTENDED | |
|------------------------------|-----------|---------------|--|
| VELOCITY HINGED KNEE PULL-UP | | | |
| XS | 13" - 15" | (33CM - 38CM) | |
| S | 15" - 17" | (38CM - 43CM) | |
| М | 16" - 18" | (41CM - 46CM) | |
| L | 18" - 20" | (46CM - 51CM) | |
| XL | 20" - 22" | (51CM - 56CM) | |
| XXL | 22" - 24" | (56CM - 61CM) | |
| XXXL | 24" - 26" | (61CM - 66CM) | |
| | | | |

| XXXL | 24" - 26" | (61CM - 66CM) |
|------------|-------------|------------------|
| VELOCITY H | HINGED KNEE | ANTERIOR CLOSURE |
| XS | 12" - 16" | (31CM - 41CM) |
| S | 14" - 18" | (36CM - 46CM) |
| М | 16" - 20" | (41CM - 51CM) |
| L | 17" - 22" | (43CM - 56CM) |
| XL | 19" - 25" | (48CM - 64CM) |
| XXL | 23" - 28" | (58CM - 71CM) |
| XXXL | 25" - 29" | (63CM - 74CM) |
| | | |



- Constructed of *Kuhl*[™] perforated neoprene to maximize patient comfort and compliance
- Range-of-Motion hinges and stop kit allows for injury specific ROM settings
- Brace can be locked at full extension
- Indicated and ideally suited for treatment of ACL, PCL, MCL and LCL sprains, as well as post-operative and non-operative meniscus injuries
- Available in anterior closure and pull-up
- Removable hinges for easy brace cleaning
- Interchangeable ext/flex stop provide 10 degrees increases in control



SHIELDS II® BRACE

- Provides dynamic control of the patella reducing discomfort & enhancing recovery
- Trakaderm® U-shaped buttress allows for treatment of multiple patella indications and fine tuning the position of the patella no matter the malalignment
- One brace fits both left & right for less Inventory
- Low Profile and extremely durable
- Independent calf and thigh straps adjust for contours of the leg
- Hinges provide medial and lateral stability can be removed for step-down application (6675)





6674 SHIELDS II® BRACE HINGED SHIELDS II® BRACE 6675

DESCRIPTION

| SIZING: | MEASURE 4" | ABOVE MID-PATELLA |
|---------|-------------|-------------------|
| XS | 10" - 11.5" | (25CM - 29CM) |
| S | 12" - 13.5" | (30CM - 34CM) |
| М | 14" - 16" | (36CM - 41CM) |
| L | 16.5" - 19" | (42CM - 48CM) |
| XL | 19.5" - 22" | (50CM - 56CM) |
| XXL | 22.5" - 25" | (57CM - 64CM) |
| | | |

| | MEASURE AT MID-PATELLA | |
|-----|------------------------|---------------|
| XS | 10" - 11" | (25CM - 28CM) |
| S | 11.5" - 13" | (29CM - 33CM) |
| М | 13.5" - 15" | (34CM - 38CM) |
| L | 15" - 17" | (38CM - 43CM) |
| XL | 17" - 19" | (43CM - 48CM) |
| XXL | 19" - 21" | (48CM - 53CM) |
| | | |



VIEW APPLICATION VIDEO T HELY-WEBER.COM

CLARENCE L. SHIELDS, M.D. Cedars-Sinai® Kerlan-Jobe Institute



KNEESIO™ BRACE

- Extremely light weight "only 8 oz."
- Flexible hinge contours to patients leg
- Breathable material keeps athlete cool while stabilizing patella tracking
- Four adjustable patella straps provides customized patellar support
- Removable/adjustable buttress for additional control where needed





7675 KNEESIO™ BRACE

C L1820 APPROVED

DESCRIPTION

| SIZING: | MEASURE 6 | " ABOVE MID-PATELLA |
|---------|-----------|---------------------|
| XS | 13" - 15" | (33CM - 38CM) |
| S | 15" - 18" | (38CM - 46CM) |
| M | 18" - 20" | (46CM - 50CM) |
| L | 20" - 23" | (50CM - 58CM) |
| XL | 23" - 25" | (58CM - 64CM) |
| XXL | 25" - 28" | (64CM - 71CM) |



VIEW APPLICATION VIDEO

WENDELL BECTON, M.D. The Orthopedic Center of St. Louis



27