

TINY CAST-AWAY®

ELBOW ORTHOSIS USPN D819219



COMFORT & INCREASED EARLY FUNCTIONALITY

5612 TINY CAST-AWAY®

INDICATIONS:

- STABLE SUPRACONDYLAR FRACTURE
- STABLE EPICONDYLAR FRACTURE
- STABLE RADIAL HEAD/NECK FRACTURE
- STABLE OLECRANON FRACTURE
- OCCULT FRACTURE
- ELBOW DISLOCATION
- ELBOW STRAIN/SPRAIN

FEATURES AND BENEFITS:

- **DUAL ROD STAYS** – PROVIDE SUPERIOR RIGIDITY AND SUPPORT WHILE REMAINING LOW PROFILE
- **LIGHT-WEIGHT, LOW-PROFILE DESIGN** – LEAVES THE ADJACENT, UNAFFECTED JOINTS UNRESTRICTED AND PROMOTES COMPLIANCE IMPROVING CLINICAL OUTCOMES
- **SIMPLE APPLICATION** – PROMOTES GOOD HYGIENE/SKIN CARE AND ELIMINATES CAST COMPLICATIONS
- **VERSATILE STRAPPING SYSTEM** – ALLOWS FOR CUSTOMIZED FIT AND COMFORT
- **MEMORY FOAM OLECRANON PAD** – PROVIDES COMFORT AND SUPPORT TO PREVENT ANY PRESSURE/IRRITATION TO VULNERABLE AREAS.

DESIGNED IN
COLLABORATION WITH:



DR. BYRON IZUKA, MD
CHILDRENS ORTHOPEDICS OF HAWAII

HELY & WEBER

BRACING INNOVATIONS

ELIMINATE THE COMPLICATIONS OF A CAST AND THE ANXIETY THAT COMES WITH THE INEVITABILITY OF SAWING IT OFF. PROMOTE GOOD HYGIENE/SKIN CARE WITH THIS SIMPLE APPLICATION.

LOW PROFILE DUAL ROD STAYS
OFFER SUPERIOR RIGIDITY AND PREMIUM SUPPORT RESULTING IN INCREASED EARLY FUNCTIONALITY.



VERSATILE STRAPPING SYSTEM
ALLOWS FOR A COMFORTABLE, CUSTOMIZED FIT THAT PATIENTS AND THEIR PARENTS WILL PREFER.



3 SIZES AVAILABLE, RIGHT/LEFT AMBIDEXTROUS
REDUCING INVENTORY NEEDS



OPTIONAL ACCESSORY:
ORTHOWEAR ELBOW (5614)
* SIZES XXS THROUGH XL



5612 ORDERING INFORMATION:

SIZE	ARM MEASUREMENT
S	6.0-7.5 (15CM-19CM)
M	7.5-8.5 (19CM-22CM)
L	8.5-10.0 (22CM-25CM)

MEASURE FROM OLECRANON TO ULNAR STYLOID

CALL TODAY! **(800) 654-3241**
EMAIL: **ORDER@HELY-WEBER.COM**

© 2021 HELY & WEBER

HELY & WEBER
BRACING INNOVATIONS