



## 220 TINY POST-OP™ KNEE

### FEATURES AND BENEFITS:

- **STRAPS:**  
Plush high-resilience foam strapping is gentle on sensitive skin.
- **HINGE:**  
Compact, durable design with easy flexion, extension, and drop-lock adjustments, plus a center lock that is less intuitive for children.
- **UPRIGHTS:**  
Composite materials provide excellent strength while contouring for a close, comfortable, and supportive fit.
- **ANTI-SLIP FEATURES:**  
Silicone-lined pads at key points prevent brace migration.
- **LIGHTWEIGHT:**  
Weighs only 1.6 lbs, ideal for kids of all sizes.
- **ADJUSTABLE:**
  - **LENGTH:** Adjusts from 14.5" to 19" with a shorter thigh upright for a natural fit.
  - **CIRCUMFERENCE:** Alligator clips allow for a simple, customized fit.
  - **STRAP TRAVELERS:** Easily repositioned by pushing and sliding.



### INDICATIONS:

- **LIGAMENTOUS INJURIES**
- **FRACTURE CARE**
- **CARTILAGINOUS INJURIES**



**LCODE  
SUGGESTED**  
L1832/L1833



### 220TINY ORDERING INFORMATION:

SIZE

PEDIATRIC UNIVERSAL FITS 11"-25"  
MEASURED 6" ABOVE MID-PATELLA

PHONE: (800) 654-3241 FAX: (800) 559-5975  
755 REGENT BLVD. PO BOX 612965 DALLAS, TX 75261-2965  
EMAIL: [ORDER@HELY-WEBER.COM](mailto:ORDER@HELY-WEBER.COM)