

# PEDIATRIC

**ORDER No. DESCRIPTION**

450 TINY TINY TITAN® WRIST  
 450-XS TITAN WRIST™ XS  
 USPN 6,893,410; 6,960,176  
 LCODS: **18308** SUGGESTED

**SIZING:**

TINY	< 4"	(10CM)
XS	4"-6"	(10CM - 15CM)

SPECIFY RIGHT OR LEFT



## TINY TITAN® WRIST

- Pediatric splint constructed of soft, comfortable orthopedic felt laminated to brush nylon exterior
- Anatomically contoured palmar stay is malleable for customization and removable for step-down rehabilitation
- Stockinette tongue eases application
- Anatomically designed for a better fit

**ORDER No. DESCRIPTION**

455 TINY TINY TITAN® THUMB  
 455-XS TITAN THUMB™ XS  
 LCODS: **18307/18307** SUGGESTED  
 USPN 6,893,410; 6,960,176  
 7,033,331; 7,056,298

**SIZING:**

TINY	< 4"	(10CM)
XS	4"-6"	(10CM - 15CM)

SPECIFY RIGHT OR LEFT



## TINY TITAN® THUMB

- Pediatric splint constructed of soft, comfortable orthopedic felt laminated to brush nylon exterior
- Stockinette tongue eases application
- Anatomically designed for a better fit

**ORDER No. DESCRIPTION**

5818-PED KUHL™ MODABBER™ WRIST  
 LONG PEDIATRIC  
 5819-PED KUHL™ MODABBER™ WRIST  
 SHORT PEDIATRIC  
 LCODS: **18308** SUGGESTED

**SIZING:**

PEDIATRIC	4" - 6"	(10CM - 15CM)
-----------	---------	---------------

SPECIFY RIGHT OR LEFT



## MODABBER™ WRIST ORTHOSIS

- Patented patient assist single strap for easy application
- Kuhl™ perforated neoprene construction provides comfortable compression
- 4 way stretch neoprene ensures brace conforms to contour of hand
- Aluminum palmar stay is removable and malleable for a custom fit

## TKO® THE KNUCKLE ORTHOSIS

**ORDER No. DESCRIPTION**

3848-XS TKO® (THE KNUCKLE ORTHOSIS)  
 3848 TKO® (THE KNUCKLE ORTHOSIS)  
 LCODS: **18307** SUGGESTED  
 USPN 7,278,980; 7,442,177;  
 7,455,650; 7,402,149;  
 7,276,039

**SIZING:**

	MEASURE CIRCUMFERENCE AT PALMAR CREASE
X-SMALL	<7" (18CM)
REGULAR	7"-10.5" (18CM-25CM)

SPECIFY RIGHT OR LEFT



- Provides positioning for fractures or injuries to metacarpals, phalanges, MCP and IP joints
- Four applications in one brace
- Straps are adjustable for fit and comfort (pain & swelling)
- Distal buddy strapping eliminates rotation of the finger
- Padded stays are malleable and allow for customized fit
- Reversible - one splint acts as a Boxer's Fracture splint (ulnar gutter), Radial Gutter splint, and Middle & Ring Finger splint


## TINY CAST-AWAY®

**ORDER No. DESCRIPTION**

5612 TINY CAST-AWAY®  
 LCODS: **18980/18982/18762** SUGGESTED

**SIZING:**

	MEASURE FROM OLECRANON TO ULNAR STYLOID
SMALL	6.0 - 7.5 (15CM - 19CM)
MEDIUM	7.5 - 8.5 (19CM - 22CM)
LARGE	8.5 - 10.0 (22CM - 25CM)



- Dual rod stays- provide superior rigidity and support while remaining low profile
- Memory foam olecranon pad - provides comfort and support to prevent any pressure or irritation to vulnerable areas
- Versatile strapping system - allows for customized fit and comfort
- Available in 3 sizes, right/left ambidextrous reducing inventory needs
- Lightweight, low profile design- promotes compliance improving clinical outcomes
- Simple application- promotes good hygiene/skin care and eliminates cast complications and cast saw anxiety



Bryon Izuka, M.D.  
 Children's Orthopaedics  
 of Hawaii LLC

