



Warnings and Instructions: Review carefully, proper application is required


⚠ Warning: This device will not prevent or eliminate risk of injury. Do Not Overtighten. If swelling, pain, skin irritation, or an unusual reaction occurs, discontinue use immediately and consult your medical professional.
 Care: Hand wash using mild soap. Rinse thoroughly. Air dry only. Do not tumble dry.




1. Align and slide thumb into thumb brace.




2. Pull tightly and secure assist strap across back of hand.



3. Stretch strap over assist strap, bring around to the palmar side of the hand, through the webspace of the thumb, and secure to hook on back of the thumb.



4. Shown with strap secured.




5. Contour the moldable aluminum buttress to support the CMC Joint comfortably.




Instructions Video

Warnings and Instructions: Review carefully, proper application is required


⚠ Warning: This device will not prevent or eliminate risk of injury. Do Not Overtighten. If swelling, pain, skin irritation, or an unusual reaction occurs, discontinue use immediately and consult your medical professional.
 Care: Hand wash using mild soap. Rinse thoroughly. Air dry only. Do not tumble dry.




1. Align and slide thumb into thumb brace.




2. Pull tightly and secure assist strap across back of hand.



3. Stretch strap over assist strap, bring around to the palmar side of the hand, through the webspace of the thumb, and secure to hook on back of the thumb.



4. Shown with strap secured.



5. Contour the moldable aluminum buttress to support the CMC Joint comfortably.



Instructions Video