

ORDER No. DESCRIPTION PHOMFIT™ WRIST HAND THUMB ORTHOSIS

447 PHOMFIT™ WRIST HAND THUMB ORTHOSIS  
 PDAC L3809 APPROVED

**SIZING:** WRIST CIRCUMFERENCE

XS	4.5" - 5.5"	(11CM - 14CM)
S	5" - 6.5"	(13CM - 17CM)
M	6" - 7.5"	(15CM - 19CM)
L	7" - 8.5"	(18CM - 22CM)
XL	8" - 9.5"	(20CM - 24CM)

SPECIFY RIGHT OR LEFT



- Aluminum removable palmar stay and radial stay allow for custom adjustment
- Stockinette tongue and dorsal loop & lock closures offer simple application
- Constructed of polyester outer shell with a tricot covered foam lining for comfort and durability
- Embedded strap guides improve fit and lower profile

ORDER No. DESCRIPTION CMC CONTROLLER™ WITH PRESSURE DOT

3801 CMC CONTROLLER™ WITH PRESSURE DOT  
 USPN D605776; 7,887,497


**SIZING:**

XS/S	4.5" - 6.5"	(11CM - 17CM)
M/L	6.5" - 10"	(17CM - 25CM)

X-LARGE SIZE ALSO AVAILABLE



- Places the thumb in the optimum position by simultaneously abducting the thumb and seating the CMC joint
- Low profile
- Universal; fits both the right and left
- Pressure dot provides support to the CMC Joint



**David M. Auerbach, M.D.**  
 Southern California Orthopedic Institute

ORDER No. DESCRIPTION SANTA BARBARA THUMB SPLINT®

3842 SANTA BARBARA THUMB SPLINT®  
 PDAC L3923/L3924 APPROVED

**SIZING:** THUMB HOLE WRIST CIRC.

UNIVERSAL	2.5" - 3.4"	4.5" - 9.5"
	(6CM - 9CM)	(11CM - 24CM)



- Constructed of Kuhl™ perforated neoprene for cool, comfortable compression
- Radial, palmar and dorsal pockets contain removable and interchangeable malleable aluminum and spring steel stays for customized support
- Wrap around design for one-handed application
- Thumb opening cut with added elastics offers close comfortable fit without constriction
- Universal. Fits right and left

**Michael Behrman, M.D.**  
 Associated Hand Surgeons

ORDER No. DESCRIPTION MODABBER™ THUMB ORTHOSIS

- Constructed of Kuhl™ perforated neoprene for cool, comfortable compression
- Brace fits left or right hand to reduce inventory
- Aluminum stay and steel spiral stay are interchangeable and removable
- Patient assist strap for easy application



5803 KUHL™ MODABBER™ THUMB ORTHOSIS  
 PDAC L3924 APPROVED

**SIZING:**

XS		
THUMB HOLE	2.1" - 2.4"	(5CM - 6CM)
WRIST CIRC.	4.5" - 5.5"	(11CM - 14CM)
REG		
THUMB HOLE	2.4" - 3"	(6CM - 8CM)
WRIST CIRC.	5" - 9.5"	(13CM - 24CM)
XL		
THUMB HOLE	3" - 3.6"	(8CM - 9CM)
WRIST CIRC.	5" - 9.5"	(13CM - 24CM)

**M. Ramin Modabber, M.D.**  
 Cedars-Sinai®  
 Kerlan-Jobe Institute

ORDER No. DESCRIPTION DEQ® THUMB ORTHOSIS

- Open stay reduces pressure over the first dorsal compartment preventing pain
- Trimmable straps and bendable stay customized fit to the patient
- Universal size ambidextrous reduced inventory needs
- Low profile
- Functional
- Comfortable



4813 DEQ® THUMB ORTHOSIS  
 PDAC L3924 APPROVED  
 PATENT PENDING

**SIZING:** THUMB HOLE WRIST CIR

UNIVERSAL	2.5" - 3.4"	4.5" - 9.5"
	(6CM - 9CM)	(11CM - 24CM)

**Steven Shin, M.D.**  
 Cedars-Sinai®  
 Kerlan-Jobe Institute

**Carey Goble, PA-C**  
 Cedars-Sinai®  
 Kerlan-Jobe Institute

**Roy Meals, M.D.**