

GELrod™ ELBOW



OPTIMAL COMFORT & PROTECTION

3837 GELrod™ ELBOW

INDICATIONS:

- CUBITAL TUNNEL SYNDROME
- OLECRANON BURSITIS/BURSECTOMY

FEATURES AND BENEFITS:

- **SPECIALLY FORMULATED GEL PAD** TO PROVIDE PROTECTION AND COMFORT, WHILE REMAINING LOW PROFILE AND COMPRESSIVE.
- **REMOVABLE ROD-STAYS** EASILY CONVERT BRACE FROM NIGHT TIME IMMOBILIZATION TO DAYTIME PROTECTION.
- **JUST TWO SIZES**, FITS RIGHT OR LEFT, KEEPS INVENTORY NEEDS TO A MINIMUM.



ROD-STAYS
REMOVED

DESIGNED IN
COLLABORATION WITH:



JAMES SAUCEDO, MD
HOUSTON METHODIST



3837 ORDERING INFORMATION:

SIZING CHART

S/M 11.0" (28CM) TO 13.0" (33CM)
L/XL 13.0" (33CM) TO 15.0" (38CM)

MEASURE WIDEST PART OF BICEP.



CALL TODAY! (800) 654-3241
EMAIL: ORDER@HELY-WEBER.COM

    © 2021 HELY & WEBER

HELY & WEBER
BRACING INNOVATIONS

Cubital Comfort™ BRACE



A NIGHTTIME COMFORT SOLUTION

3847 CUBITAL COMFORT™ BRACE

INDICATIONS:

- CUBITAL TUNNEL SYNDROME

FEATURES AND BENEFITS:

- **TWO ALUMINUM MALLEABLE STAYS**
CONTOUR TO PLACE YOUR PATIENT IN DESIRED POSITION OF ELBOW FLEXION.
- **ADJUSTABLE** TO FIT A WIDE CIRCUMFERENCE OF ARMS.
- **GLYPHIC LOCATORS** TO ASSIST IN PATIENT REAPPLICATION.
- **UNIVERSAL URETHANE STRETCHABLE STRAPS** DESIGNED FOR EASY PATIENT APPLICATION AND ADJUST FROM EITHER END OF THE STRAP TO ALLOW FOR AMBIDEXTROUS APPLICATIONS.



DESIGNED IN
COLLABORATION WITH:



DAVID M. AUERBACH, MD
SOUTHERN CALIFORNIA
ORTHOPEDIC INSTITUTE



3847 ORDERING INFORMATION:

SIZING CHART

SHORT LESS THAN 11.5" (29CM)
LONG GREATER THAN 11.5" (29CM)

MEASURE LENGTH FROM MID FOREARM TO MID HUMERUS.



CALL TODAY! (800) 654-3241
EMAIL: ORDER@HEL-Y-WEBER.COM



© 2021 HELY & WEBER

HEL-Y & WEBER

BRACING INNOVATIONS