

HELY & WEBER

BRACING INNOVATIONS

Webly® Achilles TT™

FUNCTIONAL ACHILLES PAIN RELIEF



340 WEBLY® ACHILLES TT™

A NOVEL AND INNOVATIVE BRACE THAT PROVIDES **INSTANT RELIEF** FOR ACHILLES TENDONITIS PAIN.



FEATURES AND BENEFITS:

- PATENTED POSTERIOR ENTRY GAUNTLET ENABLES SIMPLE APPLICATION
- THE POSTERIOR PANEL CONTAINS TWO SIDE BUTTRESSES THAT REST MEDIAL AND LATERAL TO THE ACHILLES TENDON FOR COMFORT
- AN ELASTIC INTERIOR BACK PANEL SEAM STRETCHES TO ENSURE ANATOMICAL POSITIONING OF THE BUTTRESSES
- THE POSTERIOR TENSION STRAP TIGHTENS, APPROXIMATING THE TWO BUTTRESSES ON BOTH SIDES OF THE ACHILLES TENDON RELIEVING PAIN
- ANTERIOR ANKLE GAUNTLET IS TIGHTENED VIA THE BOA® FIT SYSTEM FOR MICRO-ADJUSTABLE, PRECISION FIT
- THREE SIZES FIT BOTH RIGHT AND LEFT, MINIMIZING INVENTORY NEEDS

INDICATIONS:

- ACHILLES TENDONITIS
- SEVERS DISEASE



BOA® Technology integrates micro-adjustable dials, super-strong lightweight laces, and low friction lace guides onto the brace to deliver a fit solution purpose-built for performance.

AN EASY APPLICATION FOR **PAIN REDUCTION.**



The Webyl[®] Achilles TT[™] is a patented and patent pending brace that relieves pressure from irritated tissues, provides immediate comfort, reduces pain, and improves function. It accomplishes this by anatomically supporting the Achilles tendon, allowing it to stretch instead of compressing it in a way that alters its anatomy.



OPEN THE BRACE AND LOOSEN THE TENSIONING STRAP AND LACES.



SLIDE ONTO FOOT, APPROXIMATE BUTTRESSES TO ACHILLES.



TIGHTEN THE BOA[®] FIT SYSTEM TO SECURE THE FRONT.



TIGHTEN THE TENSION STRAP FOR SNUG FIT OF BUTTRESSES AGAINST ACHILLES.

DESIGNED IN COLLABORATION WITH:



JIM JABLONKA, PT

ORDERING INFORMATION (340 XS-S, 340 M-L, 340 XL-2XL):

SIZING	FITTING BY SHOE SIZE						
MEN	NARROW	REGULAR	WIDE	WOMEN	NARROW	REGULAR	WIDE
XS - S	8.5	6.5 - 8	7	XS - S	10.5	7.5 - 10	7.5
M - L	13.5	8.5 - 13	10	M - L	14.5	10.5 - 14	11
XL - 2XL	14.5	13.5 - 14	13				

CALL TODAY! (800) 654-3241
EMAIL: ORDER@HELY-WEBER.COM



© 2026 HELY & WEBER

