

HELY & WEBER®

BRACING INNOVATIONS

TKO® 2.0

USPN 12138190

TWELVE FUNCTIONS, ONE ORTHOSIS



**PATENTED VERSATILITY:
THE ULTIMATE WRIST AND
HAND ORTHOSIS FOR
UNIVERSAL APPLICATION**

4850 TKO® 2.0

FEATURES AND BENEFITS:

- **TOTAL UNIVERSALITY:**
 - **Patented Design:** Seamlessly fits both left- and right-hand applications.
 - **Adaptable Configuration:** Effortlessly transitions from a Wrist-Hand-Finger Orthosis (WHFO) to a Hand-Finger Orthosis (HFO), catering to diverse needs with a single adjustment to remove the wrist portion of the brace.
 - **Universal Sizing:** Expertly designed to accommodate most hand sizes with adjustable features, ensuring a perfect fit.
- **MULTIPLE APPLICATIONS:**
 - Both right and left hand configurations for **ulnar gutter**, **buddy splinting** (3rd/4th digits) and **radial gutter** for all-encompassing care.
- **CUSTOMIZABLE COMFORT:**
 - **Padded Malleable Stays:** Features both volar and dorsal stays that can be precisely bent to achieve optimal patient-specific positioning.
 - **Distal Buddy Strapping:** Prevents finger rotation, facilitating correct alignment and promoting effective healing.
- **EFFICIENT INVENTORY MANAGEMENT:**
 - **SKU Reduction:** Streamlines inventory with one SKU capable of covering all application needs, significantly reducing inventory costs and complexity.





ULNAR GUTTER

▪ Provides critical immobilization for the ulnar side of the wrist and hand, securing the ring and small fingers in an 'intrinsic plus' position and promoting optimal alignment and recovery.

Ideal for: Boxer's fracture, 4th/5th metacarpal (MC) fractures, metacarpophalangeal (MCP) joint injuries, proximal phalanx fractures, and ulnar nerve injuries



ULNAR GUTTER



BUDDY SPLINT

▪ Ensures stable immobilization of the ring and middle fingers, maintaining them in 'safe position'. This brace provides tailored support for less common central hand injuries.

Ideal for: 3rd/4th MC fractures, MCP joint injuries, and proximal phalanx fractures

BUDDY SPLINT



RADIAL GUTTER

▪ Offers robust protection and immobilization for the radial side, particularly the index and middle fingers, while allowing unhindered functionality of the thumb. This application maintains the fingers in an optimal position of healing.

Ideal for: 2nd/3rd MC fractures, MCP joint injuries, and proximal phalanx fractures

RADIAL GUTTER



4850 ORDERING INFORMATION:

SIZING:
REGULAR 6.0" - 11.0" (FITS MOST HANDS)
MEASURE CIRCUMFERENCE AT PALMAR CREASE

PHONE: (800) 654-3241 FAX: (800) 559-5975
755 REGENT BLVD. PO BOX 612965 DALLAS, TX 75261-2965
EMAIL: ORDER@HELY-WEBER.COM