

HELY & WEBER

BRACING INNOVATIONS

CMC Controller™ Pro

USPN D1073962

A PRACTICAL PARTNER IN CMC CARE

1001 CMC CONTROLLER™ PRO

A RIGID, OFF-THE-SHELF THUMB ORTHOSIS DESIGNED TO SUPPORT CMC ARTHRITIS WITH STABILITY AND SIMPLICITY.



FEATURES AND BENEFITS:

- **RIGID PLASTIC FRAME WITH SOFT OVERMOLDING**
Provides the perfect rigidity for supported pinching and gripping activities while maintaining durability over time.
- **SOFTER WEBSPACE CONSTRUCTION**
Offers relief in the webspace to improve comfort and compliance.
- **EMBEDDED MALLEABLE METAL STAY**
Allows for customization to enhance fit and optimize joint alignment based on patient needs.
- **NON-STRETCH STRAPS WITH NON-SLIP BACKING**
Provides secure long-lasting durability without loosening during wear.
- **THREE-POINT STRAP PANEL DESIGN**
Distributes pressure evenly across all straps to reduce pressure points and enhance overall comfort and control.

INDICATIONS:

- CMC ARTHRITIS PAIN AND INSTABILITY
- POST-OP CARE

1001 ORDERING INFORMATION (SPECIFY LEFT OR RIGHT):

SIZING	MEASURE CIRCUMFERENCE AT PALMAR CREASE	
S	6.9" - 7.8"	(18CM - 20CM)
M	7.8" - 8.9"	(20CM - 23CM)
L	8.9" - 10.6"	(23CM - 27CM)



CALL TODAY! (800) 654-3241
EMAIL: ORDER@HELY-WEBER.COM



© 2025 HELY & WEBER

