

The FingerFix™

EFFICIENCY MEETS EFFICACY



901/902 THE FINGERFIX™

OFFERS SOLUTIONS FOR MOST FINGER INJURIES WITH A SINGLE BRACE.

FEATURES AND BENEFITS:

- FULLY WELDED CONSTRUCTION REDUCES BULK AND IMPROVES COMFORT AND FUNCTIONALITY
- NON-SLIP INTERFACE MATERIAL PREVENTS MIGRATION AND MAINTAINS POSITIONING
- TRIMMABLE TO ALLOW FOR CUSTOMIZATION
- ALUMINUM STAY CAN BE BENT FOR CUSTOMIZED POSITIONING
- CAN BE APPLIED VOLARLY OR DORSALLY DEPENDING ON NEED
- 4 SIZES TO ACCOMMODATE A VARIETY OF APPLICATIONS






INDICATIONS:

- MALLET FINGER
- BOUTONNIERE DEFORMITY
- FINGER INJURY
- CRUSH INJURY
- TRIGGER FINGER



Application Video

DIP Immobilization | PIP Immobilization | DIP + PIP Immobilization | Thumb Immobilization | Trigger Finger Application



DESIGNED IN COLLABORATION WITH:



DAVID AUERBACH, MD

901/902 ORDERING INFORMATION:

901 FINGERFIX™ 4-PACK (ALL 1 SIZE)
902 FINGERFIX™ SIZE SAMPLER (ALL 4 SIZES INCLUDED)

SIZING	LENGTH OF RIGID STAY
S	1.4" (3.5cm)
M	1.8" (4.5cm)
L	2.2" (5.5cm)
XL	2.6" (6.5cm)

CALL TODAY! (800) 654-3241
EMAIL: ORDER@HELY-WEBER.COM

