

HELY & WEBER

BRACING INNOVATIONS

Liga-Loc™

LOCK IN HEALING. KEEP LIFE MOVING.

1020 LIGA-LOC™

LIGA-LOC™ DELIVERS CUSTOM-MOLDED THUMB STABILIZATION FOR UCL INJURIES, FRACTURES, AND CMC DYSFUNCTION - LOCKING IN HEALING WHILE KEEPING PATIENTS FUNCTIONAL.




FEATURES AND BENEFITS:

- **CUSTOM-MOLDED FIT:** Pressure-molded directly to the patient for intimate, secure support.
- **TARGETED IMMOBILIZATION:** Two embedded aluminum stays provide precise positioning and stabilization of the thumb.
- **FREEDOM OF MOVEMENT:** Hand-based design protects the thumb while leaving the wrist and thumb IP joint free.
- **LOW-PROFILE COMFORT:** Fully welded construction creates a smooth, sleek, and irritation-free fit.
- **SECURE HOLD:** Dual closure straps at the hand and thumb ensure consistent immobilization during activities.
- **STREAMLINED SIZING:** Two size options (S/M & L/XL) fit most patients, simplifying inventory management.

INDICATIONS:

- UCL INJURIES OF THE THUMB
- STABLE THUMB FRACTURES
- CMC ARTHRITIS OR DYSFUNCTION



EMBEDDED ALUMINUM STAYS



1020 ORDERING INFORMATION:

SIZING	MEASURE CIRCUMFERENCE AT PALMAR CREASE:
S/M	6.5" to 8.5"
L/XL	8.0" to 10.5"
	SPECIFY RIGHT OR LEFT



CALL TODAY! (800) 654-3241
EMAIL: ORDER@HELY-WEBER.COM

