

# MODABBER™ WRIST ORTHOSIS



YOUR **NEXT** PATIENT'S **FAVORITE** WRIST **BRACE**

## 5818/5819 KUHL™ MODABBER™ WRIST

### INDICATIONS:

- EXCELLENT FOR CARPAL TUNNEL SYNDROME, WRIST PAIN, STRAIN, DORSAL GANGLION CYST, AND OA OF THE WRIST

### FEATURES AND BENEFITS:

- CONSTRUCTED OF DURABLE KUHL™ PERFORATED NEOPRENE, ALLOWING MOISTURE TO ESCAPE, WHILE RETAINING HEAT AND PROVIDING COMPRESSION
- PLUSH EXTERIOR MATERIAL MAKES FOR A SLEEK APPEARANCE AND IDEAL SLEEPING BRACE
- PATIENT-ASSIST STRAP MAKES INDEPENDENT APPLICATION SIMPLE
- ALUMINUM PALMAL STAY IS REMOVABLE AND MALLEABLE FOR A CUSTOM FIT
- UNIVERSAL SIZING REDUCES INVENTORY NEEDS

DESIGNED IN  
COLLABORATION WITH:



RAMIN MODABBER, MD



### 5818 (LONG)/5819 (SHORT) ORDERING INFORMATION:

SIZE	WRIST CIRCUMFERENCE
PED	4"-6" (10-15cm)
STANDARD	4.5"-9.5" (11-24cm)
SPECIFY RIGHT OR LEFT	

CALL TODAY! (800) 654-3241  
EMAIL: [ORDER@HELY-WEBER.COM](mailto:ORDER@HELY-WEBER.COM)



© 2019 HELY & WEBER



**HELY & WEBER**  
BRACING INNOVATIONS