

HELY & WEBER

BRACING INNOVATIONS

TFO[®] The Fracture Orthosis

THE QUICK AND SIMPLE FRACTURE FIX

475 TFO[®] THE FRACTURE ORTHOSIS

A PREFABRICATED BRACE SPECIFICALLY DESIGNED AS AN ALTERNATIVE TO PLASTER AND SYNTHETIC SPLINTING MATERIALS FOR THE TREATMENT OF ACUTE HAND AND WRIST INJURIES AND FRACTURES.



USPN: 12,318,323

FEATURES AND BENEFITS:

- **THREE SIZES, FITS THE RIGHT OR LEFT HAND:**
Keeps inventory needs to a minimum.
- **EMBEDDED MALLEABLE ALUMINUM STAYS:**
Provides excellent stabilization and customizable positioning.
- **WIDE, OPEN-KNIT ELASTIC STRAPS:**
Allow for even, breathable compression and comfort.
- **MODULAR DESIGN:**
Customizable to immobilize desired joints, also excellent for step-down application.

INDICATIONS:

- ACUTE WRIST, HAND, FINGER INJURIES
- STABLE FRACTURE OF THE WRIST, HAND AND/OR FINGERS
- CRUSH INJURIES



475 ORDERING INFORMATION:

SIZING	WRIST CIRCUMFERENCE	LENGTH OF BRACE
S	5.5"-7.5" (14cm to 19cm)	9" (23cm)
M	6.5"-9.5" (17cm to 24cm)	9" (23cm)
L	7.5"-11.5" (19cm to 29cm)	9" (23cm)



CALL TODAY! (800) 654-3241
EMAIL: ORDER@HELY-WEBER.COM





WIDE, BREATHABLE,
OPEN-KNIT ELASTIC
STRAPS PROVIDE
EVEN,
COMFORTABLE
COMPRESSION TO
HELP CONTROL
EDEMA



THREE SIZES FIT
THE RIGHT *OR*
LEFT HAND,
REQUIRING ONLY
3 SKUs



MODULAR DESIGN
ALLOWS FOR
CUSTOMIZATION
TO IMMOBILIZE
DESIRED JOINTS
AND/OR FOR
STEP-DOWN
APPLICATION



A QUICKLY APPLIED WRIST/HAND IMMOBILIZER THAT TAKES YOUR PATIENT FROM ACUTE INJURY TO RETURN TO PLAY AS THEY HEAL.

THE BEST AND MOST
EFFECTIVE

HELY & WEBER
BRACING INNOVATIONS