

HELY & WEBER

BRACING INNOVATIONS

TFO™ The Fracture Orthosis

THE QUICK AND SIMPLE FRACTURE FIX

475 TFO™ THE FRACTURE ORTHOSIS

A PREFABRICATED BRACE SPECIFICALLY DESIGNED AS AN ALTERNATIVE TO PLASTER AND SYNTHETIC SPLINTING MATERIALS FOR THE TREATMENT OF ACUTE HAND AND WRIST INJURIES AND FRACTURES.



FEATURES AND BENEFITS:

- **THREE SIZES, FITS THE RIGHT OR LEFT HAND:**
Keeps inventory needs to a minimum.
- **EMBEDDED MALLEABLE ALUMINUM STAYS:**
Provides excellent stabilization and customizable positioning.
- **WIDE, OPEN-KNIT ELASTIC STRAPS:**
Allow for even, breathable compression and comfort.
- **MODULAR DESIGN:**
Customizable to immobilize desired joints, also excellent for step-down application.

INDICATIONS:

- ACUTE WRIST, HAND, FINGER INJURIES
- STABLE FRACTURE OF THE WRIST, HAND AND/OR FINGERS
- CRUSH INJURIES

475 ORDERING INFORMATION:

SIZING	WRIST CIRCUMFERENCE (MIN TO MAX)
475-S	4.5" (11cm) to 7.5" (19cm)
475-M	6.5" (17cm) to 9.5" (24cm)
475-L	7.5" (19cm) to 11.5" (29cm)



CALL TODAY! (800) 654-3241
EMAIL: ORDER@HELY-WEBER.COM



© 2025 HELY & WEBER





WIDE, BREATHABLE,
OPEN-KNIT ELASTIC
STRAPS PROVIDE
EVEN,
COMFORTABLE
COMPRESSION TO
HELP CONTROL
EDEMA



THREE SIZES FIT
THE RIGHT *OR*
LEFT HAND,
REQUIRING ONLY
3 SKUs



MODULAR DESIGN
ALLOWS FOR
CUSTOMIZATION
TO IMMOBILIZE
DESIRED JOINTS
AND/OR FOR
STEP-DOWN
APPLICATION



A QUICKLY APPLIED WRIST/HAND IMMOBILIZER THAT TAKES YOUR PATIENT FROM ACUTE INJURY TO RETURN TO PLAY AS THEY HEAL.

THE BEST AND MOST
EFFECTIVE

HELY & WEBER
BRACING INNOVATIONS