

HELY & WEBER

BRACING INNOVATIONS

CAST-AWAY®

ELBOW ORTHOSIS

USPN: 9,782,285



5613 CAST-AWAY®

INDICATIONS:

- STABLE ELBOW FRACTURES
- ELBOW DISLOCATION
- DISTAL BICEPS REPAIR AND OTHER SURGICAL INTERVENTIONS
- ELBOW INJURY

FEATURES AND BENEFITS:

- **DUAL ROD STAYS** – PROVIDE SUPERIOR RIGIDITY AND SUPPORT WHILE REMAINING LOW PROFILE
- **LIGHT-WEIGHT, LOW-PROFILE DESIGN** – LEAVES THE ADJACENT, UNAFFECTED JOINTS UNRESTRICTED AND PROMOTES COMPLIANCE IMPROVING CLINICAL OUTCOMES
- **SIMPLE APPLICATION** – PROMOTES GOOD HYGIENE/SKIN CARE AND ELIMINATES CAST COMPLICATIONS
- **VERSATILE STRAPPING SYSTEM** – ALLOWS FOR CUSTOMIZED FIT AND COMFORT
- **TELESCOPING HAND CUFF** – SUPPORTS THE WRIST AND CONTROLS FOREARM ROTATION WHEN NEEDED
- **REMOVABLE SLING STRAP** – ELIMINATES THE EARLY NEED FOR A SLING IN CONJUNCTION WITH THE BRACE
- **ADJUSTABLE HEIGHT HUMERAL CUFF** – ALLOWS FOR A CUSTOMIZED FIT

DESIGNED IN
COLLABORATION WITH:



AARON DALUISKI, MD
HOSPITAL FOR SPECIAL SURGERY, NEW YORK, NY

THE GO-TO SOLUTION FOR ELBOW IMMOBILIZATION

HELY & WEBER

BRACING INNOVATIONS



REMOVABLE SLING STRAP
OFFERS ADDITIONAL SUPPORT TO KEEP THE CAST-AWAY® COMFORTABLY STABLE AND IN A PASSIVE RESTING POSITION DURING THE RECOVERY PROCESS.



1

REMOVABLE, RIGID IMMOBILIZATION AND COMPRESSION THAT ALLOWS FOR EARLY CONTROLLED MOTION TO PREVENT RESIDUAL STIFFNESS ASSOCIATED WITH CASTING.

OPTIONAL TELESCOPING HAND CUFF
SUPPORTS THE WRIST AND CONTROLS FOREARM ROTATION WHEN NEEDED.



2

PROXIMAL AND DISTAL CUFFS ARE REMOVABLE FOR STEP-DOWN SUPPORT AS HEALING PROGRESSES AND REHABILITATION ADVANCES.



3



OPTIONAL ACCESSORY:
ORTHOWEAR ELBOW (5614)
* SIZES XXS THROUGH XL



USPN 9,782,285

5613 ORDERING INFORMATION:

SIZE	ARM CIRCUMFERENCE (IN)	H (IN)
S	8-10	8
M	10-12	9
L	12-14	10

MEASURE WIDEST PART OF FOREARM
HEIGHT OF HUMERUS INCLUDED FOR REFERENCE



CALL TODAY! **(800) 654-3241**
EMAIL: **ORDER@HELY-WEBER.COM**

