

TKO[®]

THE KNUCKLE ORTHOSIS

USPN 7,276,039



ONE SPLINT FOUR APPLICATIONS



3848 TKO[®]

INDICATIONS:

- FRACTURES OR INJURIES TO
 - METACARPALS
 - PHALANGES
 - MCP JOINTS
 - PIP JOINTS

FEATURES AND BENEFITS:

- STRAPS ADJUSTABLE FOR FIT AND COMFORT (PAIN AND SWELLING)
- DISTAL BUDDY STRAPPING ELIMINATES ROTATION OF THE FINGER
- PADDED STAYS ARE MALLEABLE FOR CUSTOM FIT TO ACHIEVE APPROPRIATE POSITIONING
- ONE SPLINT OFFERS EXCELLENT VERSATILITY BY ACTING AS A BOXER'S FRACTURE SPLINT (ULNAR GUTTER), RADIAL GUTTER SPLINT, AND MIDDLE/RING FINGER SPLINT.

DESIGNED IN COLLABORATION WITH:



DAVID GARELICK, MD
MIDWEST ORTHOPAEDICS
ILLINOIS BONE & JOINT
CHICAGO, ILLINOIS



EDWARD COTTON, OTC
COOK COUNTY HOSPITAL
CHICAGO, ILLINOIS

TKO®



**FOUR APPLICATIONS IN
ONE VERSATILE BRACE**




**BOXER'S
FRACTURE**



**3RD & 4TH
BUDDY SPLINT**



**RADIAL
GUTTER**



**DUPUYTREN'S
SPLINT**

3848 ORDERING INFORMATION:

SIZE	SIZING DESCRIPTION
3848-XS	FITS LESS THAN 7" (<18CM)
3848	FITS 7"-10" (18CM - 25CM)
3848-XL	FITS GREATER THAN 10" (>25CM)

MEASURE CIRCUMFERENCE AT PALMAR CREASE
SPECIFY RIGHT OR LEFT



CALL TODAY! (800) 654-3241
EMAIL: ORDER@HELY-WEBER.COM



HELY & WEBER

BRACING INNOVATIONS