

PDAC APPROVED L3923

SANTA BARBARA THUMB SPLINT

One Brace, Many Options

- Constructed of Kuhl™ perforated neoprene for cool, comfortable compression
- Radial, palmar, and dorsal stay pockets contain a malleable aluminum stay and a spiral stay, which can be mixed and matched for desired support
- Ideal for treating Gamekeeper's/Skier's Thumb, Basal Joint Arthritis, and Thumb Tendonitis
- Circumferential strap fits a wide range of thumb sizes
- Wrap around design for one-handed application
- Fits both the right and left hand



Stays are stackable, removable, and interchangeable



Developed in collaboration with Michael Behrman, M.D. of Associated Hand Surgeons, Santa Barbara, CA

ORDERING INFORMATION

ORDER No. DESCRIPTION / OPTIONS

3842 SANTA BARBARA THUMB SPLINT

ORDERS: (800) 654-3241

order@hely-weber.com

SIZING

THUMB HOLE SIZING
2.5"-3.375" (6.4 CM-8.6CM)

WRIST CIRCUMFERENCE
4.5"-9.5" (11.4CM-24.1CM)



HELY & WEBER

BRACING
INNOVATIONS

PDAC APPROVED L3923

MODABBER™

THUMB ORTHOSIS

US PATENT 6,142,966



- Constructed of Kuhl™ perforated neoprene for cool, comfortable compression
- Brace fits left or right hand to reduce inventory
- Aluminum stay and steel spiral stay are interchangeable and removable
- Ideal for treating Basal Joint Arthritis, Throttle Thumb, and Thumb Tendonitis
- Elasticized patient assist strap makes one handed application easy
- Taping above stay allows trimming for freedom of the IP joint



Elasticized patient assist strap makes one handed application easy



Pocket flap makes it possible to interchange and remove stays

ORDERING INFORMATION

ORDER No. DESCRIPTION / OPTIONS

5803 KUHL™ MODABBER™ THUMB ORTHOSIS

SIZING

SIZE	THUMB HOLE SIZING	WRIST CIRCUMFERENCE
XS:	2.125" - 2.375" (5.4CM - 6.0CM)	4.5" - 5.5" (11.4CM - 14.0CM)
REG:	2.375" - 3.0" (6.0CM - 7.6CM)	5.0" - 9.5" (12.7CM - 24.1CM)
XL:	3.0" - 3.625" (7.6CM - 9.2CM)	5.0" - 9.5" (12.7CM - 24.1CM)



Developed in collaboration with Dr. Modabber who is a Hand & Upper Extremity Specialist at Santa Monica Orthopaedic and Sports Medicine Group and Saint John's Health Center in Santa Monica, CA.

PDAC APPROVED L3923

NEOPRENE THUMB ORTHOSIS

New Prebent Design

- Universal size - one size fits all
- Elasticized straps for comfortable compression
- Neoprene construction provides an excellent fit of the hand and thumb
- Supports the MP and CMC joints while allowing full finger motion
- An excellent resting splint for Thumb Tendonitis and de Quervain's Tendonitis



ORDERING INFORMATION

ORDER No. DESCRIPTION / OPTIONS

3813 NEOPRENE THUMB ORTHOSIS

SIZING

THUMB HOLE SIZING	WRIST CIRCUMFERENCE
2.5"-3.375" (6.4 CM-8.6CM)	4.5"-9.5" (11.4CM-24.1CM)



Prebent Design makes an already simple application even easier

ORDERS:
(800) 654-3241
order@hely-weber.com