## THE TITAN WRIST TM USPN 6,893,410; 6,960,176



## **DIFFERENCE** CUSTOM FIT

### 450/452 TITAN WRIST™ LACING ORTHOSIS

### **FEATURES AND BENEFITS:**

- PERFORATED MATERIALS PROVIDE MOISTURE WICKING & BREATHABILITY
- LACING SYSTEM ADJUSTS TO A WIDE RANGE OF WRIST SIZES
- PALMAL AND DORSAL STAYS ARE MALLEABLE FOR CUSTOMIZATION AND REMOVABLE FOR "STEP-DOWN" REHABILITATION
- REPOSITIONABLE DORSAL STAY POD PROVIDES CUSTOM SIZING AND FIT
- ANATOMICALLY DESIGNED FOR A MORE INTIMATE FIT
- ADJUSTABLE WEB SPACE STRAP FOR COMFORT
- ON SMALL HANDS, THE EYELETS WILL CLAMSHELL OVER THE LACES FOR A SMOOTH OVERLAP
- LONGER LENGTH AVAILABLE WITH 452 TITAN WRIST™ AND FOREARM
- 450 TINY TITAN WRIST™ IS AVAILABLE FOR PEDIATRIC PATIENTS









#### **450/452 ORDERING INFORMATION:**

WRIST CIRCUMFERENCE

4.0"-6.0" (10-15cm) 4.5"-9.5" (11-24cm) 7.5"-11" (19-28cm)

XL SIZE AVAILABLE FOR THE 450 (SHORT) ONLY

CALL TODAY! (800) 654-3241 EMAIL: ORDER@HELY-WEBER.COM









f o in S © 2021 HELY & WEBER



PDAC

APPROVEL

# THE TITAN THUMB USPN 6,893,410; 6,960,176; 7,033,331; 7,056,298



## OFTEN BUT NEVER COPIED BUT DUPLICATED

### 455/459 TITAN THUMB™ ORTHOSIS

### **FEATURES AND BENEFITS:**

- PERFORATED MATERIALS PROVIDE MOISTURE WICKING & BREATHABILITY
- LACING SYSTEM ADJUSTS TO A WIDE RANGE OF SIZES
- MOVABLE STAY POD MAKES THIS BRACE UNIVERSAL AND REDUCES INVENTORY
- STOCKINETTE TONGUE OFFERS SIMPLE APPLICATION
- DESIGNED TO BE MOLDED AND CUSTOMIZED FOR PROPER POSITIONING
  - MALLEABLE OBLIQUE THUMB STAY ALLOWS FOR DESIRED THUMB POSITION
  - MALLEABLE VOLAR STAY ADJUSTS THE WRIST POSITION
  - MALLEABLE RADIAL STAY CONTOURS AND ADJUSTS THUMB ABDUCTION
  - MALLEABLE ULNAR STAY ADJUSTS TO WRIST CONTOUR AND BLOCKS DEVIATION
- ONCE SET-UP, ORTHOSIS IS EASILY REMOVED AND REAPPLIED BY PATIENT
- LONGER LENGTH AVAILABLE WITH 459 LONG TITAN THUMB™ ORTHOSIS
- 455 TINY TITAN THUMB™ IS AVAILABLE FOR PEDIATRIC PATIENTS







### **455/459 ORDERING INFORMATION:**

WRIST CIRCUMFERENCE

4.0"-6.0" (10-15cm) XL 7.5"-11" (19-28cm)
SPECIFY RIGHT OR LEFT

XL & XS SIZE AVAILABLE FOR THE 455 (SHORT) ONLY

CALL TODAY! (800) 654-3241 EMAIL: ORDER@HELY-WEBER.COM







f oin 🗡 © 2021 HELY & WEBER



