

SIZED WRIST & THUMB ORTHOSIS

REENGINEERED FOR WHEN SIZES ARE PREFERRED



434 SUEDE LACING WRIST ORTHOSIS (SHORT)
444 SUEDE LACING WRIST ORTHOSIS (LONG)

FEATURES AND BENEFITS:

- REENGINEERED LACE GUIDES PROVIDE A LOW-PROFILE AND SLEEK LOOK
- CONSTRUCTED OF DURABLE SUEDE LAMINATED TO MOISTURE MOVING POLY-FELT
- ANATOMICAL DESIGN PROVIDES IMPROVED FIT
- MALLEABLE AND REMOVABLE STAYS ALLOW FOR CUSTOMIZATION



436 PHOMFIT™ WRIST ORTHOSIS (SHORT)
437 PHOMFIT™ WRIST ORTHOSIS (LONG)

FEATURES AND BENEFITS:

- EMBEDDED STRAP GUIDES ELIMINATE BULKY D-RINGS FOR A SLEEK LOOK AND A BETTER LOW-PROFILE FIT
- COMFORTABLE NYLON AND FOAM LINING LAMINATED TO A DURABLE POLYESTER OUTER
- MALLEABLE AND REMOVABLE STAYS ALLOW FOR CUSTOMIZATION
- NO EXPOSED LATERAL HOOK PROVIDES A SNAG FREE ORTHOSIS



447 PHOMFIT™ WRIST HAND THUMB ORTHOSIS

FEATURES AND BENEFITS:

- EMBEDDED STRAP GUIDES ELIMINATE BULKY D-RINGS FOR A SLEEK LOOK AND A BETTER LOW-PROFILE FIT
- PALMAL AND RADIAL STAYS ARE REMOVABLE AND MALLEABLE FOR CUSTOMIZATION
- CONSTRUCTED OF COMFORTABLE NYLON AND FOAM LINING LAMINATED TO A DURABLE POLYESTER OUTER



SIZING: WRIST CIRCUMFERENCE

XS	4.5" - 5.5" (11CM - 14CM)
S	5" - 6.5" (13CM - 17CM)
M	6" - 7.5" (15CM - 19CM)
L	7" - 8.5" (18CM - 22CM)
XL	8" - 9.5" (20CM - 24CM)

PLEASE SPECIFY LEFT OR RIGHT

CALL TODAY! (800) 654-3241
EMAIL: ORDER@HELY-WEBER.COM



© 2019 HELY & WEBER



BRACING INNOVATIONS

THE UNO'S REENGINEERED UNIVERSAL ORTHOSIS



428 UNO WHO® WRIST HAND ORTHOSIS

US PATENT 6,960,176

FEATURES AND BENEFITS:

- EMBEDDED STRAP GUIDES ELIMINATE BULKY D-RINGS FOR A SLEEK LOOK AND A BETTER LOW-PROFILE FIT
- ULNAR STAY IS REPOSITIONABLE TO PROVIDE A CUSTOM SIZE AND FIT AND BLOCK ULNAR DEVIATION. THIS ADJUSTABILITY SIGNIFICANTLY REDUCES INVENTORY CARRYING NEEDS
- CONSTRUCTED OF COMFORTABLE NYLON AND FOAM LINING LAMINATED TO A DURABLE POLYESTER OUTER



429 UNO WHT® WRIST HAND THUMB ORTHOSIS

US PATENTS: 7,056,298; 6,960,176

FEATURES AND BENEFITS:

- EMBEDDED STRAP GUIDES ELIMINATE BULKY D-RINGS FOR A SLEEK LOOK AND A BETTER LOW-PROFILE FIT
- VOLAR AND RADIAL STAYS ARE REMOVABLE AND MALLEABLE FOR A CUSTOM FIT
- ULNAR STAY IS REPOSITIONABLE TO PROVIDE A CUSTOM SIZE AND FIT AND BLOCK ULNAR DEVIATION. THIS ADJUSTABILITY SIGNIFICANTLY REDUCES INVENTORY CARRYING NEEDS
- MALLEABLE OBLIQUE THUMB STAY ALLOWS DESIRED POSITIONING AND STABILITY OF THE THUMB
- CONSTRUCTED OF COMFORTABLE NYLON AND FOAM LINING LAMINATED TO A DURABLE POLYESTER OUTER
- ADJUSTABLE WEB SPACE STRAP ENSURES COMFORT



SIZING: WRIST CIRCUMFERENCE

UNIV 4.5" - 9.5" (11CM - 24CM)

PLEASE SPECIFY LEFT OR RIGHT

CALL TODAY! (800) 654-3241
EMAIL: ORDER@HELW-WEBER.COM



© 2019 HELY & WEBER



BRACING INNOVATIONS