

Matt Strap II

Counterforce • Focal Point Compression • Shock Absorption



The Next Generation of Counterforce Bracing for Upper Extremity

- Provides comfortable localized pressure to the forearm to reduce pain
- Malleable aluminum bolster allows for customizable pressure
- Low profile buttress creates a pseudo-attachment allowing tendon origin to rest and heal
- Non-slip interface keeps brace in place even during activity
- Universal, customizable (not a one size fits all) sizing provides individual fit with reduced inventory

Indications:

- Lateral / Medial Epicondylitis



Simple closure creates desired compression

HELY & WEBER

**BRACING
INNOVATIONS**

Matt Strap II

Counterforce • Focal Point Compression • Shock Absorption



- The strap and buttress combination reduces adverse stress on the patella tendon
- Shock absorbing buttress diffuses energy and reduces the load on the patella tendon
- Pain reduction allows the patient to resume activity
- Malleable aluminum bolster conforms to the contours of the tibia for improved function and comfort

Indications:

- Bursitis
- Osgood-Schlatter
- Runners and Jumpers Knee (Patella Tendinitis)



Low profile malleable aluminum bolster

The Next Generation of Counterforce Bracing for Lower Extremity

Ordering Information

ORDER No.	DESCRIPTION
6709	MATT STRAP II UNIVERSAL CIRCUMFERENCE: 7.0" - 20.0" (18CM - 51CM)

ORDERS: (800) 654-3241
order@hely-weber.com



*Matthew Shapiro, M.D.
Slocum Center for Orthopedics &
Sports Medicine Eugene, OR*

