

HELY & WEBER

BRACING INNOVATIONS

Kneesio™ Patellofemoral Brace

CUSTOM SUPPORT FOR ANTERIOR KNEE PAIN



7675 KNEESIO™ BRACE

MIMICKING THE BENEFITS OF TAPING, THE KNEESIO™ BRACE OFFERS LIGHT, TARGETED SUPPORT THAT MOVES WITH YOU—NO LIMITS, NO HASSLE.



FEATURES AND BENEFITS:

- LIGHTWEIGHT - ONLY 8 OZ.
- FLEXIBLE HINGE UPRIGHTS CONTOUR TO THE LEG AS WELL AS PROVIDE SUPPORT.
- BREATHABLE MATERIAL KEEPS ATHLETE COOL WHILE STABILIZING PATELLA TRACKING.
- ADJUSTABLE PATELLA STRAPS PROVIDE CUSTOMIZED PATELLAR SUPPORT.
- REMOVABLE/ADJUSTABLE BUTTRESS FOR ADDITIONAL CONTROL WHERE NEEDED.
- INTERNAL HINGES CONFORM TO THE LEG AND CONDYLE PADS PROVIDE UNPARALLELED COMFORT.
- 13" LENGTH FOR ADDED MEDIAL AND LATERAL SUPPORT.

INDICATIONS:

- PATELLAR TENDONITIS - "JUMPER'S KNEE"
- PATELLA MALALIGNMENT
- CHONDROMALACIA PATELLA
- PATELLA SUBLUXATION
- RUNNER'S KNEE

COMFORTABLE AND ADJUSTABLE

IMPROVE THE QUALITY OF **PATIENTS'** ACTIVITIES



The Kneesio™ is a game changer for patients with active lifestyles. Not only does it control pain and swelling, it maximizes patellar control, holds the patella in the proper position, and allows patients to participate in their activity of choice without a second thought. Avoid the high price tape market and slip into something way more comfortable!



REMOVABLE & ADJUSTABLE BUTTRESS CAN BE PLACED TO PROVIDE ADDITIONAL CONTROL.



THE SHEER CHARACTERISTICS OF THE BRACE MAKE IT FEEL INVISIBLE.



ELASTIC STRAPS ABSORB THE STRESSES EXERTED ON THE PATELLA, ALLOWING PAIN FREE ACTIVITIES.



LIGHTWEIGHT MATERIAL PROVIDES SUPPORT, VENTILATION, AND COMPRESSION FOR EDEMA CONTROL.



DESIGNED IN COLLABORATION WITH:



WENDELL BECTON, MD

7675 ORDERING INFORMATION (UNIVERSAL LEFT/RIGHT):

SIZING	MEASUREMENTS TAKEN AT 6" ABOVE MID-PATELLA				
XS	13" - 15"	(33CM - 38CM)	L	20" - 23"	(50CM - 58CM)
SM	15" - 18"	(38CM - 46CM)	XL	23" - 25"	(58CM - 64CM)
M	18" - 20"	(46CM - 50CM)	XXL	25" - 28"	(64CM - 71CM)

CALL TODAY! (800) 654-3241
EMAIL: ORDER@HELY-WEBER.COM



© 2025 HELY & WEBER

