

RAPID ZAP™

ANKLE ORTHOSIS



MAXIMUM PROTECTION **IN** ALL ZONES



318 RAPID ZAP™

INDICATIONS:

- STRAINS AND SPRAINS OF THE ANKLE
- PROPHYLACTIC USE

FEATURES AND BENEFITS:

- **UNIQUE THREE ZONE LACING SYSTEM** ENSURES UNIFORM LACE TIGHTNESS AND RAPID DONNING AND DOFFING.
- **FIGURE 8 STRAPS** MIMIC ANKLE TAPING IN A PREFERRED NEUTRAL POSITION PROVIDING RESISTANCE TO INVERSION/EVERSION ANKLE INJURIES YET NOT RESTRICTING FLEXION.
- **STRAPS ARE ANGLED** TO ENSURE "UNDER HEEL" APPLICATION FOR MAXIMUM SUPPORT AND COMFORT.
- **STURDY CONSTRUCTION AND SEAMLESS DESIGN** OFFER ADDED COMFORT WITH TOP IN CLASS QUALITY, AND DURABILITY.
- **ELASTIC CUFF CLOSURE** SECURES LACES AND STRAPS WHILE OFFERING COMPRESSION.

FIXING THE SPEED LACER ONE ZONE AT A TIME...

INDEPENDENT
ADJUSTMENT
FOR FOREFOOT
ENSURES
UNIFORM LACE
TIGHTNESS




Zone 1

MIDFOOT
STRAP
SECURES
TALUS AND
PROVIDES
GREATER
STABILITY




Zone 2

TWO POINT
ADJUSTING
TRAVELER
STRAP
CREATES A
WIDE THROAT
FOR A LARGE,
LONG, OR
SWOLLEN
FOOT/ANKLE
TO EASILY PASS
THROUGH



Zone 3



APPLICATION
VIDEO

318 ORDERING INFORMATION:

SIZING:	FITTING BY SHOE SIZE						
	NARROW	REGULAR	WIDE	WOMEN	NARROW	REGULAR	WIDE
MEN							
XXS	6.5	5.5-6	5	XXS	7.5	6.5-7	6
XS	7.5	6.5-7	6	XS	9	7.5-8.5	7
S	8.5	7.5-8	7	S	10.5	8-10	7.5
M	11.5	8.5-11	8	M	12.5	10.5-12	10
L	13.5	11.5-13	10	L	14.5	12-14	11
XL	14.5	13.5-14	12				
XXL		14.5-16	13				

CALL TODAY! (800) 654-3241
EMAIL: ORDER@HELY-WEBER.COM



© 2022 HELY & WEBER

HELY & WEBER[®]
BRACING INNOVATIONS