

HINGED **RAPIDZAP**™ MAXIMUM PROTECTION IN SECONDS



318H HINGED RAPID ZAP™

INDICATIONS

- STRAINS AND SPRAINS OF THE ANKLE
- WHEN MORE SUPPORT IS NEEDED THAN A STANDARD LACE UP BRACE

FEATURES

- STURDY BALLISTIC NYLON ASSEMBLY PROVIDES DURABILITY AND LONGEVITY
- SEAMLESS CONSTRUCTION ELIMINATES POTENTIAL IRRITATION
- LOW PROFILE HINGE PROVIDES RIGIDITY AND SUPPORT TO THE ANKLE
- CLOSE FITTING DESIGN SLIPS EASILY INTO THE SHOE



HELY & WEBER
BRACING INNOVATIONS

NEXT GENERATION OF THE ZAP FAMILY



LARGE THROAT

- OPENS FOR LONG FEET AND LARGE ANKLES. GET THE CORRECT BRACE FROM THE BEGINNING. NO NEED TO UPSIZE FOR A SWOLLEN ANKLE

HINGE ASSEMBLY

- UNIVERSAL LEFT/RIGHT HINGE AND FOOTPLATE ASSEMBLY EXTENDS FROM MIDFOOT TO HINDFOOT FOR MAXIMUM CONTROL AND COMFORT
- ACCOMODATES TO THE PATIENT'S ANATOMY

FIGURE 8 STRAPS

- MIMIC ANKLE TAPING TO SECURE THE ANKLE MORTISE AND CONTROL INVERSION/EVERSION

THREE ZONE SPEED LACER DESIGN

- NO NEED TO TIGHTEN THE LACES AT EACH EYELET
- THREE SPEED LACING TABS ALLOW FAST INDEPENDENT ADJUSTMENTS TO THE FORE-FOOT, MIDFOOT AND ANKLE



ZONE 1



ZONE 2



ZONE 3



318H HINGED RAPID ZAP™

SIZING: FITTING BY SHOE SIZE

WOMEN	5-6	6.5-8	8.5-10	10.5-12	12.5-14	
NARROW	XS	S	M	L		
REGULAR	XS	S	M	L		
WIDE	S	M	L	L	XL	
MEN	6-7	7.5-9	9.5-11	11.5-13	13.5-15	15+
NARROW	XS	S	M	L	XL	
REGULAR	XS	S	M	L	XL	XXL
WIDE	S	M	L	XL	XXL	XXL



APPLICATION
VIDEO

HELY & WEBER[®]
BRACING INNOVATIONS

CALL TODAY! (800) 654-3241
EMAIL: ORDER@HELY-WEBER.COM



© 2020 HELY & WEBER

ORDERING
INFO