

# 'LASTIC LUMBAR SUPPORT WITH INSERT

## Lumbosacral Supports



- Heavy elastic construction ensures comfort
- Elastic pocket accommodates thermoplastic insert
- Double side pulls cross in front for added support
- Supports the lower back and encourages normal range of motion

### Order No. Description / Options

102 'Lastic Lumbar Support for Moldable Insert  
106 Moldable Insert

Measure widest part of hip/waist:

**XS:** 32" - 35" (81.3cm - 88.9cm)    **Sm:** 36" - 41" (91.4cm - 104.1cm)  
**Med:** 42" - 45" (106.7cm - 114.3cm)    **Lg:** 46" - 50" (116.8 cm - 127.0 cm)  
**XL:** 51" - 55" (129.5 cm - 139.7 cm)    **XXL:** 56" - 60" (142.2cm - 152.4cm)



Thermoplastic Insert



**HELY & WEBER**  
*Innovative Supports*

1185 East Main Street, Santa Paula, California 93060  
California: 800.221.5465 National: 800.654.3241  
Fax: 800.559.5975 [www.hely-weber.com](http://www.hely-weber.com)