

THE GUS ARM SLING



THE WAY ALL SLINGS SHOULD BE DESIGNED

500 THE GUS ARM SLING

FEATURES AND BENEFITS:

- THE GUS' POSTERIOR GUSSET IS A COMFORT PANEL THAT ANATOMICALLY ALLOWS THE SLING TO SEAT THE ARM AND ELBOW
- HIGH RESILIENCE FOAM STRAP IS ADJUSTABLE AND COMFORTABLE
- DEEP POCKET DESIGN ACCOMMODATES THE LARGEST OF CASTS AND DRESSINGS
- TWO ADJUSTMENT POINTS ATTACHED TO O-RINGS ALLOW SLING TO BETTER FIT DIFFERING PATIENT BUILDS
- OPEN END WITH SOFT THUMB LOOP ADDS COMFORT AND ELIMINATES SLIPPAGE
- NO NEED TO PAINFULLY SHRUG YOUR SHOULDER OR FLEX YOUR ELBOW TO ENTER SLING



500 ORDERING INFORMATION:

SIZE	ARM MEASUREMENT (LENGTH)
XS	8"-10" (20CM-25CM)
S	10"-12" (25CM-30CM)
M	12"-14" (30CM-36CM)
L	14"-16" (36CM-41CM)
XL	16"-18" (41CM-46CM)

MEASURE FROM OLECRANON TO KNUCKLES
7" HEIGHT FOR ALL SIZES

CALL TODAY! (800) 654-3241

EMAIL: ORDER@HELY-WEBER.COM



© 2021 HELY & WEBER

HELY & WEBER

BRACING INNOVATIONS

THE GUS SHOULDER IMMOBILIZER



UNPARALLELED PATIENT COMFORT & SUPPORT


501 THE GUS SHOULDER IMMOBILIZER

FEATURES AND BENEFITS:

- THE GUS' POSTERIOR GUSSET IS A COMFORT PANEL THAT ANATOMICALLY ALLOWS THE SLING TO SEAT THE ARM AND ELBOW
- HIGH RESILIENCE FOAM STRAPS ARE ADJUSTABLE AND COMFORTABLE
- DEEP POCKET DESIGN ACCOMMODATES THE LARGEST OF CASTS AND DRESSINGS
- TWO ADJUSTMENT POINTS ATTACHED TO O-RINGS FOR EACH STRAP ALLOWING SLING TO BETTER FIT DIFFERING PATIENT BUILDS
- OPEN END WITH SOFT THUMB LOOP ADDS COMFORT AND ELIMINATES SLIPPAGE



POSTERIOR GUSSET AND
IMMOBILIZATION STRAP ON O-RING
TRAVELER PROVIDES OPTIMAL FIT



THUMB LOOP REST AND
WRIST STRAP TRAVELER

501 ORDERING INFORMATION:

SIZE	ARM MEASUREMENT (LENGTH)
XS	8"-10" (20CM-25CM)
S	10"-12" (25CM-30CM)
M	12"-14" (30CM-36CM)
L	14"-16" (36CM-41CM)
XL	16"-18" (41CM-46CM)

MEASURE FROM OLECRANON TO KNUCKLES
7" HEIGHT FOR ALL SIZES

CALL TODAY! (800) 654-3241

EMAIL: ORDER@HELY-WEBER.COM



© 2021 HELY & WEBER



HELY & WEBER

BRACING INNOVATIONS