

SHOULDER SOLUTIONS

SHOULDER IMMOBILIZERS



PDAC L3670 APPROVED

THE GUS SHOULDER IMMOBILIZER (501)

A TRUSTED CHOICE FOR SECURE AND COMFORTABLE SHOULDER SUPPORT

FEATURES AND BENEFITS:

- The posterior gusset allows the arm to anatomically sit within the sling
- Padded shoulder strap is built from high resiliency foam offering maximum comfort
- O-Ring adjustment points allow sling to better fit a variety of girths
- Traveler strap ensures best waist strap fit
- Straps include alligator hook for easy strap length adjustment
- Sizing available: XS, S, M, L, XL



PDAC L3670 APPROVED

THE GUS SEMI-UNIVERSAL (501)

ALL THE TRUSTED SUPPORT, LESS INVENTORY

FEATURES AND BENEFITS:

- All the same features and benefits of **the GUS Shoulder Immobilizer** with an adjustable sling body that provides semi-universal sizing reducing inventory needs to just two sizes (S/M & L/XL).



PDAC L3670 APPROVED

VAPOR IMMOBILIZATION SLING (504)

A PREMIUM SHOULDER SUPPORT PROVIDING LONG-TERM COMFORT AND BREATHABILITY

FEATURES AND BENEFITS:

- Constructed of NoSwet® material, combined with a soft moisture wicking liner provides cool comfort
- Padded shoulder strap is built from high resiliency foam offering maximum comfort
- Quick release buckles allow for easy donning/doffing
- The posterior gusset allows the arm to anatomically sit within the sling
- Straps include alligator hook for easy strap length adjustment
- Sizing available: S, M, L, XL

SHOULDER SOLUTIONS

ABDUCTION SLINGS



PDAC L3670 APPROVED

SHOULDER HOULDER II™ (514B)

EXACTLY WHAT YOU NEED WHEN REHABILITATION REQUIRES IMMOBILIZATION IN ABDUCTION

FEATURES AND BENEFITS:

- Constructed of NoSwet® material, combined with a soft moisture wicking liner provides cool comfort
- Padded shoulder strap is built from high resiliency foam offering maximum comfort
- Quick release buckles allow for easy donning/doffing
- The posterior gusset allows the arm to anatomically sit within the sling
- Straps include alligator hook for easy strap length adjustment
- Sizing available: S, M, L, XL



PDAC L3670 APPROVED

SHOULDER HOULDER II™ UNIVERSAL (514B-U)

RELIABLE ABDUCTION SUPPORT, REDUCED INVENTORY NEEDS

FEATURES AND BENEFITS:

- All the same features and benefits of the **Shoulder Houlder II™** with an adjustable sling body that provides universal sizing with just one SKU.



SIZING AND ORDERING INFORMATION

THE GUS (501)

SIZING:	MEASURE FROM OLECRANON TO KNUCKLES
XS	8" - 10"
S	10" - 12"
M	12" - 14"
L	14" - 16"
XL	16" - 18"

THE GUS SEMI-UNIVERSAL (501)

SIZING:	MEASURE FROM OLECRANON TO KNUCKLES
S/M	10" - 14"
L/XL	14" - 18"

VAPOR (504)

SIZING:	MEASURE FROM OLECRANON TO KNUCKLES
S	11.5" - 13"
M	13.5" - 14.5"
L	15" - 16"
XL	16.5" - 17.5"

SHOULDER HOULDER II™ (514B)

SIZING:	MEASURE FROM OLECRANON TO KNUCKLES
S	11.5" - 13"
M	13.5" - 14.5"
L	15" - 16"
XL	16.5" - 17.5"

SHOULDER HOULDER II™ UNIVERSAL (514B-U)

SIZING:	MEASURE FROM OLECRANON TO KNUCKLES
UNIVERSAL	11.5" - 17.5"

PHONE: (800) 654-3241 FAX: (800) 559-5975
 755 REGENT BLVD. PO BOX 612965 DALLAS, TX 75261-2965
 EMAIL: ORDER@HELY-WEBER.COM