

# The Hand Cuff Fracture Orthosis


1010

**Warnings and Instructions: Review carefully, proper application is required**


**⚠ Warning:** This device will not prevent or eliminate risk of injury. Do Not Overtighten. If swelling, pain, skin irritation, or an unusual reaction occurs, discontinue use immediately and consult your medical professional.

This device should not be worn by persons with known allergies to neoprene.


Care: Hand wash using mild soap. Rinse thoroughly. Air dry only. Do not tumble dry.




1. Align edge of brace with thenar crease and distal palmer crease as shown.



2. Bend brace around hand to form a good fit.



3. Pull and attach straps to brace at back of hand. Note: For left handed application straps must be twisted.



4. Continue to form brace for improved fit. Completed application shown.

# The Hand Cuff Fracture Orthosis


1010

**Warnings and Instructions: Review carefully, proper application is required**


**⚠ Warning:** This device will not prevent or eliminate risk of injury. Do Not Overtighten. If swelling, pain, skin irritation, or an unusual reaction occurs, discontinue use immediately and consult your medical professional.

This device should not be worn by persons with known allergies to neoprene.


Care: Hand wash using mild soap. Rinse thoroughly. Air dry only. Do not tumble dry.




1. Align edge of brace with thenar crease and distal palmer crease as shown.



2. Bend brace around hand to form a good fit.



3. Pull and attach straps to brace at back of hand. Note: For left handed application straps must be twisted.



4. Continue to form brace for improved fit. Completed application shown.