



Warnings and Instructions: Review carefully, proper application is required

⚠ Warning: This device will not prevent or eliminate risk of injury. Do Not Overtighten. If swelling, pain, skin irritation, or an unusual reaction occurs, discontinue use immediately and consult your medical professional. This device should not be worn by persons with known allergies to neoprene.


Care: Hand wash using mild soap. Rinse thoroughly. Air dry only. Do not tumble dry.



1. Insert thumb into brace.



2. Mold the primary stay to the contour of the palm and around the back of the hand.



3. Firmly affix the closure straps so that the brace is snug.



4. Bend/mold the oblique thumb stay to position the thumb as desired.



5. Firmly affix the thumb closure straps.



6. To remove the brace, undo the hand strap, slide hand/thumb out of the brace. To re-apply simply slide hand/thumb back in the brace and affix hand strap.



© 2025 Weber Orthopedics Inc.
DBA Hely & Weber all rights reserved



Application Video




Rev A 8/25


Warnings and Instructions: Review carefully, proper application is required

⚠ Warning: This device will not prevent or eliminate risk of injury. Do Not Overtighten. If swelling, pain, skin irritation, or an unusual reaction occurs, discontinue use immediately and consult your medical professional. This device should not be worn by persons with known allergies to neoprene.


Care: Hand wash using mild soap. Rinse thoroughly. Air dry only. Do not tumble dry.



1. Insert thumb into brace.



2. Mold the primary stay to the contour of the palm and around the back of the hand.



3. Firmly affix the closure straps so that the brace is snug.



4. Bend/mold the oblique thumb stay to position the thumb as desired.



5. Firmly affix the thumb closure straps.



6. To remove the brace, undo the hand strap, slide hand/thumb out of the brace. To re-apply simply slide hand/thumb back in the brace and affix hand strap.



© 2025 Weber Orthopedics Inc.
DBA Hely & Weber all rights reserved



Application Video



Rev A 8/25