

## **TKO Fitting Quick Guide:**

### **Boxers Fracture:**

Left splint fits left hand 4th and 5th metacarpals

Right splint fits right hand 4th and 5th metacarpals

### **Splinted Middle/Ring Fingers:**

Left splint fits left hand 3rd and 4th metacarpals

Right splint fits right hand 3rd and 4th metacarpals

### **Radial Gutter:**

Left splint fits right hand 2nd and 3rd metacarpals

Right splint fits left hand 2nd and 3rd metacarpals