



Warnings and Instructions: Review carefully, proper application is required

⚠ Warning: This device will not prevent or eliminate risk of injury. Do Not Overtighten. If swelling, pain, skin irritation, or an unusual reaction occurs, discontinue use immediately and consult your medical professional.


Care: Hand wash using mild soap. Rinse thoroughly. Air dry only. Do not tumble dry. For single patient use only.




1. Position hand in brace.




2. Adjust the stay pod as needed into a comfortable position on the side of the wrist.



3. Pull all three straps through guide slots.



4. Pull and fasten straps.



5. Trim excess strap material.

U.S. Patent No. 6960176

HELLY & WEBER


© 2018 Weber Orthopedics Inc.
DBA Helly & Weber all rights reserved

1185 E Main St., Santa Paula, California 93060
California: 800-221-5465, National: 800-654-3241
International: 805-525-4244, fax: 805-933-2348
U.S. fax: 800-559-5975, www.helly-weber.com


Warnings and Instructions: Review carefully, proper application is required

⚠ Warning: This device will not prevent or eliminate risk of injury. Do Not Overtighten. If swelling, pain, skin irritation, or an unusual reaction occurs, discontinue use immediately and consult your medical professional.


Care: Hand wash using mild soap. Rinse thoroughly. Air dry only. Do not tumble dry. For single patient use only.




1. Position hand in brace.




2. Adjust the stay pod as needed into a comfortable position on the side of the wrist.



3. Pull all three straps through guide slots.



4. Pull and fasten straps.



5. Trim excess strap material.

U.S. Patent No. 6960176

HELLY & WEBER

© 2018 Weber Orthopedics Inc.
DBA Helly & Weber all rights reserved

1185 E Main St., Santa Paula, California 93060
California: 800-221-5465, National: 800-654-3241
International: 805-525-4244, fax: 805-933-2348
U.S. fax: 800-559-5975, www.helly-weber.com