



Warnings and Instructions: Review carefully, proper application is required


⚠ Warning: This device will not prevent or eliminate risk of injury. Do Not Overtighten. If swelling, pain, skin irritation, or an unusual reaction occurs, discontinue use immediately and consult your medical professional.
Care: Hand wash using mild soap. Rinse thoroughly. Air dry only. Do not tumble dry.




1. Position hand in brace.




2. Position the palm stay below the midhand crease, and the knuckles above the edge of the brace.




3. Detach if necessary and align the lacing stay pod with the middle knuckle.



4. Secure the thumb strap.



5. Pull and secure the lacing strap.




6. Completed application.

U.S. Patent No: 6893410, 6960176




Warnings and Instructions: Review carefully, proper application is required


⚠ Warning: This device will not prevent or eliminate risk of injury. Do Not Overtighten. If swelling, pain, skin irritation, or an unusual reaction occurs, discontinue use immediately and consult your medical professional.
Care: Hand wash using mild soap. Rinse thoroughly. Air dry only. Do not tumble dry.




1. Position hand in brace.




2. Position the palm stay below the midhand crease, and the knuckles above the edge of the brace.




3. Detach if necessary and align the lacing stay pod with the middle knuckle.



4. Secure the thumb strap.



5. Pull and secure the lacing strap.



6. Completed application.

U.S. Patent No: 6893410, 6960176

