

Ulnar Abutment Wrist Brace


5820

Warnings and Instructions: Review carefully, proper application is required

⚠ Warning: This device will not prevent or eliminate risk of injury. Do Not Overtighten. If swelling, pain, skin irritation, or an unusual reaction occurs, discontinue use immediately and consult your medical professional.


This device should not be worn by persons with known allergies to neoprene.

Care: Hand wash using mild soap. Rinse thoroughly. Air dry only. Do not tumble dry.




1. Place thumb in thumb hole, line up stay along palm of wrist, and adjust gel pad insert against the Ulnar Styloid.


Center piece may be removed.



2. Pull up and fasten dorsum flap on top of gel assist strap.



3. Pull and fasten long strap around wrist, applying pressure as needed.



4. Completed application.



Scan to watch application video

Ulnar Abutment Wrist Brace


5820

Warnings and Instructions: Review carefully, proper application is required

⚠ Warning: This device will not prevent or eliminate risk of injury. Do Not Overtighten. If swelling, pain, skin irritation, or an unusual reaction occurs, discontinue use immediately and consult your medical professional.


This device should not be worn by persons with known allergies to neoprene.

Care: Hand wash using mild soap. Rinse thoroughly. Air dry only. Do not tumble dry.




1. Place thumb in thumb hole, line up stay along palm of wrist, and adjust gel pad insert against the Ulnar Styloid.


Center piece may be removed.



2. Pull up and fasten dorsum flap on top of gel assist strap.



3. Pull and fasten long strap around wrist, applying pressure as needed.



4. Completed application.



Scan to watch application video