


Pediatric Shoulder Immobilizing Sling


501PED

Warnings and Instructions: Review carefully, proper application is required


⚠ Warning: This device will not prevent or eliminate risk of injury. Do Not Overtighten. If swelling, pain, skin irritation, or an unusual reaction occurs, discontinue use immediately and consult your medical professional.
 Care: Hand wash using mild soap. Rinse thoroughly. Air dry only. Do not tumble dry.




1. Place arm in sling pouch.



2. Pull the shoulder strap over the shoulder, through the double rings and secure.



3. Locate the waist strap and bring to front.




4. Pull waist strap through ring and secure.

Pediatric Shoulder Immobilizing Sling


501PED

Warnings and Instructions: Review carefully, proper application is required


⚠ Warning: This device will not prevent or eliminate risk of injury. Do Not Overtighten. If swelling, pain, skin irritation, or an unusual reaction occurs, discontinue use immediately and consult your medical professional.
 Care: Hand wash using mild soap. Rinse thoroughly. Air dry only. Do not tumble dry.




1. Place arm in sling pouch.



2. Pull the shoulder strap over the shoulder, through the double rings and secure.



3. Locate the waist strap and bring to front.



4. Pull waist strap through ring and secure.

HELW&WEBER

© 2017 Weber Orthopedics Inc.
 DBA Hely & Weber all rights reserved

1185 E Main St., Santa Paula, California 93060
 California: 800-221-5465, National: 800-654-3241
 International: 805-525-4244, fax: 805-933-2348
 U.S. fax: 800-559-5975, www.hely-weber.com

CE
EC REP

MDSS GmbH, Schiffgraben 41
 30175 Hannover, Germany

Rev A 3/15

HELW&WEBER

© 2017 Weber Orthopedics Inc.
 DBA Hely & Weber all rights reserved

1185 E Main St., Santa Paula, California 93060
 California: 800-221-5465, National: 800-654-3241
 International: 805-525-4244, fax: 805-933-2348
 U.S. fax: 800-559-5975, www.hely-weber.com

CE
EC REP

MDSS GmbH, Schiffgraben 41
 30175 Hannover, Germany

Rev A 3/15