



Sling & Swathe (Universal)

523



Warnings and Instructions: Review carefully, proper application is required

⚠ Warning: This device will not prevent or eliminate risk of injury. Do Not Overtighten. If swelling, pain, skin irritation, or an unusual reaction occurs, discontinue use immediately and consult your medical professional.
Care: Hand wash using mild soap. Rinse thoroughly. Air dry only. Do not tumble dry.

The sling & swathe immobilizes the shoulder and can be used on either shoulder. Separate each component:



Slide arm through the sling. Place each strap over the left and right shoulders. Criss-cross the straps across the back. Pull the strap ends under the arm and fasten each strap to the front of the sling. Ensure that the straps are secured close to the elbow and wrist for proper support. Optionally, adjust the length of each straps by moving the metal buckle for a better fit.



Attach one end of the swathe to the front of the sling, just above the arm. Wrap the swathe around the outside of the injured arm and the body, under the unaffected arm. Secure the strap to the front of the swathe.




Sling & Swathe (Universal)

523



Warnings and Instructions: Review carefully, proper application is required

⚠ Warning: This device will not prevent or eliminate risk of injury. Do Not Overtighten. If swelling, pain, skin irritation, or an unusual reaction occurs, discontinue use immediately and consult your medical professional.
Care: Hand wash using mild soap. Rinse thoroughly. Air dry only. Do not tumble dry.

The sling & swathe immobilizes the shoulder and can be used on either shoulder. Separate each component:



Slide arm through the sling. Place each strap over the left and right shoulders. Criss-cross the straps across the back. Pull the strap ends under the arm and fasten each strap to the front of the sling. Ensure that the straps are secured close to the elbow and wrist for proper support. Optionally, adjust the length of each straps by moving the metal buckle for a better fit.



Attach one end of the swathe to the front of the sling, just above the arm. Wrap the swathe around the outside of the injured arm and the body, under the unaffected arm. Secure the strap to the front of the swathe.