



Warnings and Instructions: Review carefully, proper application is required


⚠ Warning: This device will not prevent or eliminate risk of injury. Do Not Overtighten. If swelling, pain, skin irritation, or an unusual reaction occurs, discontinue use immediately and consult your medical professional. This device should not be worn by persons with known allergies to neoprene. Care: Hand wash using mild soap. Rinse thoroughly. Air dry only. Do not tumble dry.




1. Loosen the outside straps and center the brace over the spine.




2. Pull and fasten the front hook and loop closure.




3. Center and attach the front panel over the abdomen.



4. Pull outside straps and fasten them to the front panel.




5. Completed application, front.




6. Completed application, back.

Warnings and Instructions: Review carefully, proper application is required


⚠ Warning: This device will not prevent or eliminate risk of injury. Do Not Overtighten. If swelling, pain, skin irritation, or an unusual reaction occurs, discontinue use immediately and consult your medical professional. This device should not be worn by persons with known allergies to neoprene. Care: Hand wash using mild soap. Rinse thoroughly. Air dry only. Do not tumble dry.




1. Loosen the outside straps and center the brace over the spine.




2. Pull and fasten the front hook and loop closure.




3. Center and attach the front panel over the abdomen.



4. Pull outside straps and fasten them to the front panel.



5. Completed application, front.



6. Completed application, back.