


# Auerbach Mallet Splint


600

## Warnings and Instructions: Review carefully, proper application is required


**⚠ Warning:** This device will not prevent or eliminate risk of injury. Do Not Overtighten. If swelling, pain, skin irritation, or an unusual reaction occurs, discontinue use immediately and consult your medical professional. This device should not be worn by persons with known allergies to neoprene. Care: Hand wash using mild soap. Rinse thoroughly. Air dry only. Do not tumble dry.




1. The DIP joint must be kept in a straight position at all times. Holding the fingertip on the edge of a counter, place the splint on the nail side of the finger with the rubber surface in contact with the skin.




2. Apply the splint with the transverse strap closer to the hand and the angled strap over the nail.



3. Secure the transverse strap.



4. Secure the angled strap. The straps should hold the splint firmly, but should not be too tight as to interfere with circulation to the fingertip.



5. Completed application. The PIP joint should move freely.

U.S. Patent No. D515216

**HELY & WEBER**

© 2014 Weber Orthopedics Inc.  
DBA Hely & Weber all rights reserved

1185 E Main St., Santa Paula, California 93060  
California: 800-221-5465, National: 800-654-3241  
International: 805-525-4244, fax: 805-933-2348  
U.S. fax: 800-559-5975, www.hely-weber.com

**CE**  
EC REP

MDSS GmbH, Schiffgraben 41  
30175 Hannover, Germany


Rev C 02/14

# Auerbach Mallet Splint


600

## Warnings and Instructions: Review carefully, proper application is required


**⚠ Warning:** This device will not prevent or eliminate risk of injury. Do Not Overtighten. If swelling, pain, skin irritation, or an unusual reaction occurs, discontinue use immediately and consult your medical professional. This device should not be worn by persons with known allergies to neoprene. Care: Hand wash using mild soap. Rinse thoroughly. Air dry only. Do not tumble dry.




1. The DIP joint must be kept in a straight position at all times. Holding the fingertip on the edge of a counter, place the splint on the nail side of the finger with the rubber surface in contact with the skin.




2. Apply the splint with the transverse strap closer to the hand and the angled strap over the nail.



3. Secure the transverse strap.



4. Secure the angled strap. The straps should hold the splint firmly, but should not be too tight as to interfere with circulation to the fingertip.



5. Completed application. The PIP joint should move freely.

U.S. Patent No. D515216

**HELY & WEBER**

© 2014 Weber Orthopedics Inc.  
DBA Hely & Weber all rights reserved

1185 E Main St., Santa Paula, California 93060  
California: 800-221-5465, National: 800-654-3241  
International: 805-525-4244, fax: 805-933-2348  
U.S. fax: 800-559-5975, www.hely-weber.com

**CE**  
EC REP

MDSS GmbH, Schiffgraben 41  
30175 Hannover, Germany

Rev C 02/14