

## CONTACT

CALIFORNIA:  
1185 East Main Street  
Santa Paula, CA 93060

TEXAS:  
755 Regent Blvd., Suite 250  
DFW Airport, TX 75261

NATIONAL:  
800.654.3241 or 800.221.5465  
Fax: 800.559.5975

INTERNATIONAL:  
805.525.4244  
Fax: 805.933.2348

[WWW.HEL-Y-WEBER.COM](http://WWW.HEL-Y-WEBER.COM)



OUTCOMES MATTER

HELY & WEBER<sup>®</sup>

HELY & WEBER<sup>®</sup>

HELY & WEBER<sup>®</sup>

COMPANY  
PROFILE

John Hely and Jim Weber have over 96 years combined experience in the orthopedic and sports medicine industry. They started their partnership in 1992 with their main goal, "To Develop, Manufacture, and Distribute quality braces." The company evolved from a local Southern California distribution company to an international manufacturing and distribution corporation. Throughout Hely & Weber's history we have proudly been family owned and operated.

Hely & Weber's philosophy is to satisfy the clinical needs of patients with high quality, exclusive products. Our desire is to be adaptable to the ever-changing needs of the medical community; seeking out unique state-of-the-art product designs to satisfy the vast array of patient needs. Our employees are committed to our uncompromising standards for quality and customer satisfaction.

All of our products are put through rigorous quality testing before they are released to the marketplace. When you purchase our products you will find an unparalleled energy and commitment at Hely & Weber that matches your passion for patient care.

OUR MISSION

"We develop and sell high quality innovative supports and braces. Our commitment to quality, comfort, durability, fit and function sets industry standards, ensuring total customer satisfaction and confidence in our products. All of our company resources are directed, focused and linked to this purpose."

ORDERING  
INFORMATION

ORDERING POLICY:  
Hely & Weber proudly accepts orders from practitioners, medical professionals and their staff.

TO PLACE AN ORDER:  
National: (800) 654-3241 or (800) 221-5465  
International: (805) 525-4244 or Fax: (805) 933-2348  
Fax: (800) 559-5975 - faxable order form available upon request  
E-mail: order@hely-weber.com

MAILING ADDRESS:  
California: Hely & Weber, P.O. Box 832, Santa Paula, CA 93061-0832  
Texas: Hely & Weber, P.O. Box 612956, Dallas, TX 75261-2956

REIMBURSEMENT CODES  
Please refer to our website, www.hely-weber.com, for the latest reimbursement codes and PDAC approval letters.

TABLE OF  
CONTENTS

|                    |                      |    |
|--------------------|----------------------|----|
| UPPER<br>EXTREMITY | WRIST / HAND / THUMB | 4  |
|                    | HAND / FINGER        | 11 |
|                    | ELBOW                | 14 |
|                    | SHOULDER             | 17 |

|       |       |    |
|-------|-------|----|
| BACKS | SPINE | 20 |
|-------|-------|----|

|                    |                            |    |
|--------------------|----------------------------|----|
| LOWER<br>EXTREMITY | GLOBAL KNEE                | 23 |
|                    | HINGED KNEE                | 26 |
|                    | HINGED PATELLA STABILIZERS | 28 |
|                    | KNEE SLEEVES               | 30 |
|                    | IMMOBILIZERS               | 30 |
|                    | PATELLA STABILIZERS        | 31 |
|                    | PATELLA TENDONITIS STRAPS  | 32 |
|                    | ANKLE / FOOT               | 35 |

|           |                             |    |
|-----------|-----------------------------|----|
| PEDIATRIC | WRIST / HAND / THUMB        | 37 |
|           | ELBOW / SHOULDER / CLAVICAL | 38 |
|           | KNEE / ANKLE / FOOT         | 40 |

|                         |  |    |
|-------------------------|--|----|
| COMMERCIAL<br>CUSTOMERS | CUSTOM LOGO PROGRAM / CUSTOM BRACE PROGRAM | 41 |
|                         | TERMS & POLICIES                           | 42 |

UPPER  
EXTREMITY



TITAN WRIST™ LACING ORTHOSIS

- Perforated materials for breathability
- Adjustable lacing system handles a wide range of wrist sizes
- Palmar and dorsal stays are malleable for customization and removable for step-down rehabilitation
- Repositionable and malleable dorsal stay pod provides custom sizing and fit
- Anatomically designed for a better fit
- Adjustable web space strap for comfort



| ORDER No. | DESCRIPTION                     |
|-----------|---------------------------------|
| 450       | TITAN WRIST™                    |
| 452       | TITAN WRIST™ & FOREARM          |
|           | PDAC L3908 APPROVED             |
|           | USPN 6,893,410; 6,960,176       |
| SIZING:   | WRIST CIRCUMFERENCE:            |
| REG       | 4.5" - 9.5" (11.4CM - 24.1CM)   |
| XS        | 4.0" - 6.0" (10.2CM - 15.2CM)   |
| XL        | 7.5" - 11.0" (19.1CM - 27.9CM)  |
|           | XL SIZES AVAILABLE FOR 450 ONLY |
|           | SPECIFY RIGHT OR LEFT           |

SUEDE LACING WRIST ORTHOSIS

- Constructed of perforated black suede laminated to moisture moving poly-felt
- Anatomically designed for a better fit
- Lacing design provides for easy secure closures
- Aluminum palmar and dorsal stays are malleable and removable



| ORDER No. | DESCRIPTION                   |
|-----------|-------------------------------|
| 434       | LACING WRIST ORTHOSIS         |
| 444       | LACING WRIST & FOREARM        |
|           | PDAC L3908 APPROVED           |
| SIZING:   | WRIST CIRCUMFERENCE           |
| XS        | 4.5" - 5.5" (11.4CM - 14.0CM) |
| S         | 5.0" - 6.5" (12.7CM - 16.5CM) |
| M         | 6.0" - 7.5" (15.2CM - 19.1CM) |
| L         | 7.0" - 8.5" (17.8CM - 21.6CM) |
| XL        | 8.0" - 9.5" (20.3CM - 24.1CM) |
|           | SPECIFY RIGHT OR LEFT         |



ORDER No.

DESCRIPTION

5818

KUHL™ MODABBER™ WRIST LONG

5819

KUHL™ MODABBER™ WRIST SHORT

PDAC

L3908

APPROVED

SIZING:

WRIST CIRCUMFERENCE

XS

4.0" - 6.0" (10.2CM - 15.2CM)

STANDARD

4.5" - 9.5" (11.4CM - 24.1CM)

SPECIFY RIGHT OR LEFT



M. Ramin Modabber, M.D.

Cedars-Sinai®

Kerlan-Jobe Institute



5819



5818

MODABBER™ WRIST ORTHOSIS

- Kuhl™ neoprene provides compression and comfort
- Patient assist strap for easy application
- Soft plush outer material ensures nighttime comfort
- Aluminum palmar stay is removable and malleable for a custom fit
- Premium outer loop delivers long lasting durability

ORDER No.

DESCRIPTION

428

UNO WHO® WRIST HAND ORTHOSIS

PDAC

L3908

APPROVED

USPN

6,960,176

SIZING:

WRIST CIRCUMFERENCE

4.5" - 9.5" (11.4CM - 24.1CM)

SPECIFY RIGHT OR LEFT



UNO WHO® WRIST HAND ORTHOSIS

- Aluminum palmar stay is removable and malleable for a custom fit
- Movable and malleable stay pod provides custom sizing and fit
- Stockinette tongue and dorsal loop & lock closures offer simple application
- Constructed of polyester outer shell with a tricot covered foam lining for comfort and durability
- Embedded strap guides improve fit and lower profile

ORDER No.

DESCRIPTION

3846

GEL WRIST WRAP (GEL INCLUDED)

PDAC

L3908

APPROVED

SIZING:

WRIST CIRCUMFERENCE

4.5" - 9.5" (11.4CM - 24.1CM)

SPECIFY RIGHT OR LEFT



GEL WRIST WRAP

- Contoured Vitamin E infused mineral oil gel pad is an excellent shock absorber, aids in scar maturation, and eases the discomfort of pillar pain
- Patient assist strap for easy application
- Universal sizing reduces inventory
- Aluminum palmar stay is removable and malleable for a custom fit

ORDER No.

DESCRIPTION

438

UNIVERSAL WRIST ORTHOSIS

439

UNIVERSAL WRIST & FOREARM

PDAC

L3908

APPROVED

SIZING:

WRIST CIRCUMFERENCE

4.5" - 9.5" (11.4CM - 24.1CM)

SPECIFY RIGHT OR LEFT



UNIVERSAL WRIST ORTHOSIS

- Constructed of comfortable orthopedic felt laminated to brushed nylon
- Palmar and dorsal malleable aluminum stays may be molded for a custom fit or removed for stepdown rehabilitation
- Universal sizing reduces inventory

ORDER No.

DESCRIPTION

436

PHOMFIT™ WRIST ORTHOSIS SHORT

437

PHOMFIT™ WRIST ORTHOSIS LONG

PDAC L3908 APPROVED

SIZING:

WRIST CIRCUMFERENCE

XS

4.5" - 5.5" (11.4CM - 14.0CM)

S

5.0" - 6.5" (12.7CM - 16.5CM)

M

6.0" - 7.5" (15.2CM - 19.1CM)

L

7.0" - 8.5" (17.8CM - 21.6CM)

XL

8.0" - 9.5" (20.3CM - 24.1CM)

SPECIFY RIGHT OR LEFT

437

436

PHOMFIT™ WRIST ORTHOSIS

- Comfortable foam and nylon tricot lining, laminated to a sturdy polyester outer fabric, provides comfort and durability
- Aluminum palmar stay is removable and malleable for a custom fit
- Stockinette tongue and dorsal hook & loop closures offer simple application
- No exposed hook
- Embedded strap guides improve fit and lower profile

ORDER No.

DESCRIPTION

5820

ULNAR ABUTMENT BRACE

PDAC L3908 APPROVED

SIZING:

UNIVERSAL

SPECIFY RIGHT OR LEFT

VIEW APPLICATION VIDEO AT HELY-WEBER.COM

ULNAR ABUTMENT BRACE

- Treats ulnar sided wrist pain through dorsal buttressing
- Reduces pain and helps stabilize the DRUJ
- Unique gel pad offers comfort and multiple configurations

Paul Pace, M.D.  
San Antonio  
Orthopaedic  
Group

ORDER No.

DESCRIPTION

2801

"SQUEEZE" ULNAR COMPRESSION WRAP

SIZING:

UNIVERSAL

VIEW APPLICATION VIDEO AT HELY-WEBER.COM

"SQUEEZE" ULNAR COMPRESSION WRAP

- Treats ulnar sided wrist pain with focal point compression
- Provides focal point compression to avoid tourniquet effect of circumferential wraps
- Perforated to let moisture escape
- Easy application
- Universal; fits both the right and left wrist

Jamie McMillan,  
OTR, CHT

ORDER No.

DESCRIPTION

455

TITAN THUMB™ ORTHOSIS

459

LONG TITAN THUMB™ ORTHOSIS

PDAC L3807/L3809 APPROVED

USPN 6,893,410; 6,960,176  
7,033,331; 7,056,298

SIZING:

WRIST CIRCUMFERENCE

REG

4.5" - 9.5" (11CM - 24CM)

XS

4" - 6" (10CM - 15CM)

XL

7.5" - 11" (19CM - 28CM)

XS & XL SIZES AVAILABLE FOR 455 ONLY  
SPECIFY RIGHT OR LEFT

459

455

TITAN THUMB™ ORTHOSIS

- Perforated materials for breathability
- Adjustable lacing pulls easily fit a wide range of wrist sizes
- Movable stay pod makes this brace universal and reduces inventory
- Stockinette tongue offers simple application
- Designed to be molded and customized for proper positioning
  - Malleable oblique thumb stay allows desired thumb position
  - Malleable volar stay adjusts the wrist position
  - Malleable radial stay contours and adjust thumb abduction
  - Malleable ulnar stay adjusts to wrist contour and blocks deviation
- Once set-up, easily removed and reapplied by patients

ORDER No.

DESCRIPTION

429

UNO WHT® (WRIST HAND THUMB ORTHOSIS)

PDAC L3807/L3809 APPROVED

USPN: 7,056,298; 6,960,176

SIZING:

WRIST CIRCUMFERENCE

4.5" - 9.5" (11CM - 24CM)

SPECIFY RIGHT OR LEFT

UNO WHT® - WRIST HAND THUMB

- Repositionable and malleable ulnar stay pod provides custom sizing and fit
- Stockinette tongue and dorsal loop & lock closures offer simple application
- Constructed of polyester outer shell with a tricot covered foam lining for comfort and durability
- Malleable oblique thumb stay allows desired thumb position and stability
- Ulnar stay contours and blocks deviation
- Adjustable web space strap ensures comfort
- Embedded strap guides improve fit and lower profile

Jamie McMillan,  
OTR, CHT



ORDER No.

DESCRIPTION

447

PHOMFIT™ WRIST HAND THUMB ORTHOSIS

PDAC

L3809

APPROVED

SIZING:

WRIST CIRCUMFERENCE

|    |             |               |
|----|-------------|---------------|
| XS | 4.5" - 5.5" | (11CM - 14CM) |
| S  | 5" - 6.5"   | (13CM - 17CM) |
| M  | 6" - 7.5"   | (15CM - 19CM) |
| L  | 7" - 8.5"   | (18CM - 22CM) |
| XL | 8" - 9.5"   | (20CM - 24CM) |

SPECIFY RIGHT OR LEFT

PHOMFIT™ WRIST HAND THUMB ORTHOSIS

- Aluminum removable palmar stay and radial stay allow for custom adjustment
- Stockinette tongue and dorsal loop & lock closures offer simple application
- Constructed of polyester outer shell with a tricot covered foam lining for comfort and durability
- Embedded strap guides improve fit and lower profile



ORDER No.

DESCRIPTION

3801

CMC CONTROLLER™ WITH PRESSURE DOT

USPN D605776; 7,887,497

SIZING:

|      |             |               |
|------|-------------|---------------|
| XS/S | 4.5" - 6.5" | (11CM - 17CM) |
| M/L  | 6.5" - 10"  | (17CM - 25CM) |

X-LARGE SIZE ALSO AVAILABLE

CMC CONTROLLER™ WITH PRESSURE DOT

- Places the thumb in the optimum position by simultaneously abducting the thumb and seating the CMC joint
- Low profile
- Universal; fits both the right and left
- Pressure dot provides support to the CMC Joint



FLUOROSCAN FEMALE, AGE 70 CONDITION: CMC ARTHRITIS

No BRACE

BRACE WITH CMC CONTROLLER™



David M. Auerbach, M.D.

Southern California Orthopedic Institute



ORDER No.

DESCRIPTION

3842

SANTA BARBARA THUMB SPLINT®

PDAC

L3923/L3924

APPROVED

SIZING:

THUMB HOLE

WRIST CIRC.

|           |             |               |
|-----------|-------------|---------------|
| UNIVERSAL | 2.5" - 3.4" | 4.5" - 9.5"   |
|           | (6CM - 9CM) | (11CM - 24CM) |

SANTA BARBARA THUMB SPLINT®

- Constructed of Kuhl™ perforated neoprene for cool, comfortable compression
- Radial, palmar and dorsal pockets contain removable and interchangeable malleable aluminum and spring steel stays for customized support
- Wrap around design for one-handed application
- Thumb opening cut with added elastics offers close comfortable fit without constriction
- Universal. Fits right and left



Michael Behrman, M.D.

Associated Hand Surgeons



ORDER No.

DESCRIPTION

5803

KUHL™ MODABBER™ THUMB ORTHOSIS

PDAC

L3924

APPROVED

SIZING:

|             |             |               |
|-------------|-------------|---------------|
| XS          |             |               |
| THUMB HOLE  | 2.1" - 2.4" | (5CM - 6CM)   |
| WRIST CIRC. | 4.5" - 5.5" | (11CM - 14CM) |

REG

|             |           |               |
|-------------|-----------|---------------|
| THUMB HOLE  | 2.4" - 3" | (6CM - 8CM)   |
| WRIST CIRC. | 5" - 9.5" | (13CM - 24CM) |

XL

|             |           |               |
|-------------|-----------|---------------|
| THUMB HOLE  | 3" - 3.6" | (8CM - 9CM)   |
| WRIST CIRC. | 5" - 9.5" | (13CM - 24CM) |

MODABBER™ THUMB ORTHOSIS

- Constructed of Kuhl™ perforated neoprene for cool, comfortable compression
- Brace fits left or right hand to reduce inventory
- Aluminum stay and steel spiral stay are interchangeable and removable
- Patient assist strap for easy application



M. Ramin Modabber, M.D.

Cedars-Sinai® Kerlan-Jobe Institute



ORDER No.

DESCRIPTION

4813

DEQ® THUMB ORTHOSIS

PDAC

L3924

APPROVED

PATENT PENDING

SIZING:

THUMB HOLE

WRIST CIR

|           |             |               |
|-----------|-------------|---------------|
| UNIVERSAL | 2.5" - 3.4" | 4.5" - 9.5"   |
|           | (6CM - 9CM) | (11CM - 24CM) |

DEQ® THUMB ORTHOSIS

- Open stay reduces pressure over the first dorsal compartment preventing pain
- Trimable straps and bendable stay customized fit to the patient
- Universal size ambidextrous reduced inventory needs
- Low profile
- Functional
- Comfortable



Steven Shin, M.D.

Cedars-Sinai® Kerlan-Jobe Institute



Carey Goble, PA-C

Cedars-Sinai® Kerlan-Jobe Institute



Roy Meals, M.D.





ORDER No.

DESCRIPTION

2804

CMC CONTROLLER PLUS®  
PDAC L3924 APPROVED  
USPN 7,854,714

SIZING:

S/MFITS MOST FEMALES

L/XLFITS MOST MALES

SPECIFY RIGHT OR LEFT

VIEW APPLICATION VIDEO  
AT HELY-WEBER.COM

- Low profile design offers pain relief and positioning for function
- Completely Trimable™ allowing customization for each patient
- Indicator arrow for effective strap placement
- The malleable stay hugs and lifts the head of the metacarpal while simultaneously seating the basilar joint, aiding in pain relief
- Patient assist strap makes for simple one-handed application
- Two sizes right and left reduces inventory

ORDER No.

DESCRIPTION

600

AUERBACH MALLET SPLINT  
PDAC L3927 APPROVED  
PACK OF 3

SIZING:

MEASURE FROM FINGER TIP TO 2ND JOINT (PIP)

S1.4" 3.5CM

M1.8" 4.5CM

L2.2" 5.5CM

XL2.6" 6.5CM

VIEW APPLICATION VIDEO  
AT HELY-WEBER.COM

- Designed to be a mallet splint, but can be used to immobilize the PIP and/or DIP joints from a volar or dorsal approach
- Brush nylon with hook & loop contact closure
- Supports the DIP joint in a neutral, extended position to prevent mallet deformity
- Malleable aluminum stay can be bent to hyperextend the DIP joint
- Allows full motion of the PIP joint
- Soft, scaled rubber underside prevents migration and slippage, even when wet
- Easy to wash and re-apply
- Packaged in units of three

David M. Auerbach, M.D.  
Southern California  
Orthopedic Institute

ORDER No.

DESCRIPTION

4849-SHORT

SHORT XTENDER®  
4849-LONG

LONG XTENDER®  
PDAC L3923 APPROVED  
USPN D810310

SIZING:

LENGTH OF BRACE

SHORT8" (20CM) FITS MOST FEMALES

LONG8.5" (22CM) FITS MOST MALES

VIEW APPLICATION VIDEO  
AT HELY-WEBER.COM

- A post-operative/post-procedural solution to Dupuytren's splinting
- Fits both left and right hands in two sizes for minimal inventory requirement
- Malleable hand pan allows desired positioning
- Can be applied over light and heavy dressings

David Garelick, M.D.  
Illinois Bone & Joint

ORDER No.

DESCRIPTION

3848

TKO® (THE KNUCKLE ORTHOSIS)  
PDAC L3807/L3809 APPROVED

4848

TKO® PREBENT  
PDAC L3809 APPROVED  
USPN: 7,442,177; 7,455,650  
7,402,149; 7,278,980;  
7,276,039

SIZING:

MEASURE CIRCUMFERENCE AT PALMAR CREASE

X-SMALL< 7" (18CM)

REGULAR7" - 10" (18CM - 25CM)

X-LARGE> 10" (25CM)

4848 NOT AVAILABLE IN XS

SPECIFY RIGHT OR LEFT

- Provides positioning for fractures or injuries to the metacarpals, proximal phalanx, and MCP joints
- Four applications in one brace
- Straps are adjustable for fit and comfort (pain & swelling)
- Distal buddy strapping eliminates rotation of the finger
- Padded stays are malleable and allow for customized fit


David Garelick, M.D.  
Illinois Bone & Joint

Edward Cotton, OTC

ORDER No. DESCRIPTION

THE HANDCUFF®

1010




**SIZING:** MEASURE THE CIRCUMFERENCE OF HAND BELOW KNUCKLES

|   |                          |
|---|--------------------------|
| S | 6.5" - 8.5" (17CM-22CM)  |
| M | 8.5" - 9.5 (22CM-24CM)   |
| L | 9.5" - 11.0" (24CM-28CM) |



- Customizable metacarpal cuff used for stable metacarpal fractures or step-down support following cast removal
- Wrist and fingers are free to promote functional use
- Compression molded technology creates excellent conformability providing a comfortable intimate fit to protect the injury site
- Three sizes; fits right and left hand




**Dr. Steven Shin, MD**  
Cedars-Sinai®  
Kerlan-Jobe Institute



ORDER No. DESCRIPTION

THE HANDCUFF® I PLUS

1048




**SIZING:** MEASURE THE CIRCUMFERENCE OF HAND BELOW KNUCKLES

|      |                         |
|------|-------------------------|
| S/M  | 6.5" - 8.5" (17CM-22CM) |
| L/XL | 8.5" - 10" (22CM-25CM)  |

SPECIFY RIGHT OR LEFT



- Unique compression molded technology for customized fit and position
- Malleable hand cuff with adjustable dorsal stay
- Comfortable and easy to don/doff
- Optional finger extension stay and strap
- Two sizes



**James Saucedo, M.D.**  
Houston Methodist

Not Pictured:  
**Abel Cavazos,**  
OTR, CHT  
The Hand Center  
of San Antonio




MTC FRACTURE BRACE®

ORDER No. DESCRIPTION

- A fresh alternative to painstaking casting and splinting procedures
- Controls forearm pronation/supination, allows elbow flexion/extension
- Simple and fast application
- Easy to adjust and step-down
- Can be applied to arm in neutral, pronated, and supinated positions
- Comfortable to wear especially versus a cast or custom-made splint




640




**MTC FRACTURE BRACE®**

650



**MTC FRACTURE BRACE® W/THUMB SPICA**

USPN: 9,782,285



| SIZING: | WIDTH       | LENGTH    |
|---------|-------------|-----------|
| S/M     | 9.7" (25CM) | 8" (20CM) |
| L/XL    | 12" (30CM)  | 9" (23CM) |



| WRIST COMPONENT DIMENSIONS |            |                 |
|----------------------------|------------|-----------------|
|                            | 640 WRIST  | 650 THUMB SPICA |
| XS                         | 7" (18CM)  | 8" (20CM)       |
| REGULAR                    | 8" (20CM)  | 9" (23CM)       |
| LONG                       | 10" (25CM) | 11" (28CM)      |

| SIZING:    | ARM LENGTH                                     | ARM CIRC.                      | COMPONENTS USED     |
|------------|--|--------------------------------|---------------------|
|            | MEASURE FROM OLECRANON TO DISTAL PALMAR CREASE | MEASURE WIDEST PART OF FOREARM |                     |
| 640-S-8-XS | 10"-14" (25-36CM)                              | 7.5"-11" (19-28CM)             | 8"XS WRIST : S CUFF |
| 640-S-8    | 10"-15" (25-38CM)                              | 9"-11" (23-28CM)               | 8" WRIST : S CUFF   |
| 640-S-10   | 11"-17" (28-43CM)                              | 9"-11" (23-28CM)               | 10" WRIST : S CUFF  |
| 640-L-8    | 11"-16" (28-41CM)                              | 11"-14" (28-36CM)              | 8" WRIST : L CUFF   |
| 640-L-10   | 11"-18" (28-46CM)                              | 11"-16" (28-41CM)              | 10" WRIST : L CUFF  |
| 650-S-8-XS | 11"-15" (28-38CM)                              | 7.5"-11" (19-28CM)             | 8"XS THUMB : S CUFF |
| 650-S-8    | 12"-16" (30-41CM)                              | 9"-11" (23-28CM)               | 8" THUMB : S CUFF   |
| 650-S-10   | 12"-18" (30-46CM)                              | 9"-11" (23-28CM)               | 10" THUMB : S CUFF  |
| 650-L-8    | 12"-17" (30-43CM)                              | 11"-14" (28-36CM)              | 8" THUMB : L CUFF   |
| 650-L-10   | 13"-18" (33-46CM)                              | 11"-16" (28-41CM)              | 10" THUMB : L CUFF  |

SPECIFY RIGHT OR LEFT

**Jamie McMillan,**  
OTR, CHT





ORDER No. DESCRIPTION

CAST-AWAY®

5613 CAST-AWAY®  
PDAC L3980, 3982 APPROVED  
USPN 9,782,285

5614 ORTHOWEAR ELBOW  
SIZES XXS THRU XL  
PDAC L3995 APPROVED

| CAST-AWAY™<br>SIZING: | MEASURE<br>WIDEST PART<br>OF FOREARM | HEIGHT OF<br>HUMERUS INCL.<br>FOR REFERENCE |
|-----------------------|--------------------------------------|---|
|                       | ARM CIR.                             | HEIGHT                                      |
| S                     | 8" - 10"                             | 8"  |
| M                     | 10" - 12"                            | 9"  |
| L                     | 12" - 14"                            | 10"   |

| ORTHOWEAR<br>ELBOW<br>SIZING: | MEASURE CIRCUMFERENCE AT<br>WIDEST PART OF FOREARM |
|-------------------------------|--|
|-------------------------------|--|

|     |                             |
|-----|-----------------------------|
| XXS | 7.5" - 8.5" (19CM - 22CM)   |
| XS  | 8.5" - 10" (22CM - 25CM)    |
| S   | 9.5" - 11" (24CM - 28CM)    |
| M   | 11" - 12.5" (28CM - 32 CM)  |
| L   | 11.5" - 13.5" (29CM - 34CM) |
| XL  | 12" - 14" (30CM - 36CM)     |



VIEW APPLICATION VIDEO  
AT HELY-WEBER.COM



Aaron Daluiski, M.D.  
Hospital for Special  
Surgery, New York, NY



- Dual Rod Stays provide superior rigidity and support while remaining low profile
- Lightweight, low profile design leaves the adjacent, unaffected joints unrestricted and promotes compliance improving clinical outcomes
- Simple Application promotes good hygiene/skin care and eliminates cast complications
- Versatile Strapping Systems allows for customized fit and comfort
- Telescoping Hand Cuff supports the wrist and controls forearm rotation when needed
- Removable Sling Strap eliminates the early need for a sling in conjunction with the brace
- Adjustable Height Humeral Cuff allows for a customized fit
- Removable Sling Strap offers additional support to keep the Cast-Away® comfortably stable and in a passive resting position during the recovery process
- Removable, rigid immobilization and compression that allows for early controlled motion to prevent residual stiffness associated with casting
- Proximal and Distal Cuffs are removable for step-down support as healing progresses and rehabilitation advances
- Fits right or left arm



5614  
OPTIONAL ACCESSORY: ORTHOWEAR  
ELBOW SIZES XXS THROUGH XL

CUBITAL COMFORT™ BRACE

ORDER No. DESCRIPTION

- Ideal for treatment of Cubital Tunnel Syndrome
- Two aluminum malleable stays contour to place your patient in desired position of elbow flexion
- Adjustable to fit a wide circumference of arms
- Glyphic locators to assist in simple patient reapplication
- Fits right or left arm



3847 CUBITAL COMFORT™ BRACE  
PDAC L3762 APPROVED

| SIZING: | MEASURE LENGTH FROM MID<br>FOREARM TO MID HUMERUS |
|---------|---|
| SHORT   | <11 1/2" (29CM)                                   |
| LONG    | >11 1/2" (29CM)                                   |



VIEW APPLICATION VIDEO  
AT HELY-WEBER.COM

David M. Auerbach, M.D.  
Southern California  
Orthopedic Institute



MATT STRAP™ II UNIVERSAL

ORDER No. DESCRIPTION

- Provides comfortable localized pressure to the forearm to reduce pain
- Malleable aluminum bolster allows for customizable pressure
- Low profile buttress creates a pseudo-attachment allowing tendon origin to rest and heal
- Universal, customizable (not a one size fits all) sizing provides individual fit with reduced inventory



6709 MATT STRAP™ II UNIVERSAL

| SIZING: | ENCIRCLE THE AREA OF<br>APPLICATION TO MEASURE<br>FITS CIRCUMFERENCE OF 7"- 20" (18CM - 51CM) |
|---------|---|
|---------|---|

Matthew Shapiro, M.D.  
Slocum Center for  
Orthopedics & Sports  
Medicine



TENNIS ELBOW - WITH PRESSURE PAD

ORDER No. DESCRIPTION

- Low profile design
- Neoprene with 2-sided nylon outer covering
- Pressure pad & elastic strap provide compression without constriction



3838 TENNIS ELBOW -  
PRESSURE PAD

| SIZING: | MEASURE THE WIDEST PART<br>OF FOREARM |
|---------|---------------------------------------|
| XS      | 8" - 9" (20CM - 23CM)                 |
| S       | 9" - 10" (23CM - 25CM)                |
| M       | 10" - 11" (25CM - 28CM)               |
| L       | 11" - 12" (28CM - 31CM)               |
| XL      | 12" - 14" (31CM - 36CM)               |
| XXL     | 14" - 16" (36CM - 41CM)               |

ORDER No.

DESCRIPTION

501

THE GUS SHOULDER IMMOBILIZER SLING

PDAC L3670 APPROVED

500

GUS ARM SLING

PDAC L3670 APPROVED

SIZING:

MEASURE FROM OLECRANON TO KNUCKLES 7" HEIGHT

|    |           |               |
|----|-----------|---------------|
| XS | 8" - 10"  | (20CM - 25CM) |
| S  | 10" - 12" | (25CM - 30CM) |
| M  | 12" - 14" | (30CM - 36CM) |
| L  | 14" - 16" | (36CM - 41CM) |
| XL | 16" - 18" | (41CM - 46CM) |





501

500

- The GUS' posterior GUSSET is an intimate comfort panel that anatomically allows the sling to seat the arm and elbow
- Padded pressure sensitive straps are adjustable and comfortable
- Deep pocket design accommodates the largest of casts and dressings
- Two adjustment points attached to O-rings for each strap allow sling to better fit differing patient's girth
- Open end with soft thumb loop adds comfort and eliminates slippage
- Waist strap attaches at two slider straps that seek their best attachment location
- Straps include alligator hook for easy strap length adjustment

ORDER No.

DESCRIPTION

504

VAPOR IMMOBILIZATION SLING

PDAC L3670 APPROVED

SIZING:

MEASURE FROM OLECRANON TO KNUCKLES

|    |               |               |
|----|---------------|---------------|
| S  | 11.5" - 13"   | (29CM - 33CM) |
| M  | 13.5" - 14.5" | (33CM - 37CM) |
| L  | 15" - 16"     | (37CM - 41CM) |
| XL | 16.5" - 17.5" | (42CM - 45CM) |



- Constructed of our NoSwet® material, combined with a soft moisture quick-wicking liner to provide your patient with premium comfort and support
- Body and shoulder straps provide quick adjustment, easy application, and padded comfort
- Quick release buckles allow for easy application and removal
- A posterior gusset comfortably holds the elbow in the desired position without providing unnecessary pressure
- Posterior O-ring on sling accommodates patient allowing for better fit and comfort
- Straps include alligator hook for easy strap length adjustment

17

ORDER No.

DESCRIPTION

514

SHOULDER HOULDER II™

PDAC L3670 APPROVED

524

SHOULDER HOULDER II™ NEUTRAL

PDAC L3670 APPROVED

SIZING:

MEASURE FROM OLECRANON TO KNUCKLES

|    |               |               |
|----|---------------|---------------|
| S  | 11.5" - 13"   | (29CM - 33CM) |
| M  | 13.5" - 14.5" | (33CM - 37CM) |
| L  | 15" - 16"     | (37CM - 41CM) |
| XL | 16.5" - 17.5" | (42CM - 45CM) |



514

524

- Constructed of our NoSwet® material, combined with a soft moisture quick-wicking liner to provide your patient with comfort and support
- 15° and 30° Abduction Pillows allow for ambidextrous and custom patient fitting
- Body and shoulder strap provide quick adjustment, easy application, and padded comfort
- Quick release buckles allow for easy application and removal
- The posterior gusset comfortably cradles the elbow in the desired position without providing unnecessary pressure
- Straps include alligator hook for easy strap length adjustment

KNAPP SAK® II - SHOULDER ABDUCTION KIT

ORDER No.

DESCRIPTION

510

KNAPP SAK® II

PDAC L3670 APPROVED

USPN 7,189,213; 7,300,410

SIZING:

UNIVERSAL



Thomas P. Knapp, M.D.

Cedars-Sinai®

Kerlan-Jobe Institute



- Maintains shoulder in neutral position allowing approximately 30 degrees of shoulder abduction
- Prevents excessive internal rotation of the shoulder
- Color coded for easy application

ORDER No.

DESCRIPTION

661

CLAVICLE SPLINT

PDAC L3650 APPROVED

SIZING:

MEASURE FROM CENTER OF SHOULDER BLADES AROUND SHOULDER

|    |     |     |      |       |
|----|-----|-----|------|-------|
|    | MIN | MAX | MIN  | MAX   |
| XS | 15" | 21" | 38CM | 53CM  |
| S  | 19" | 28" | 48CM | 71CM  |
| M  | 24" | 34" | 61CM | 86CM  |
| L  | 30" | 38" | 76CM | 97CM  |
| XL | 34" | 40" | 86CM | 102CM |



- Comfortable foam construction provides ideal immobilization
- Contoured urethane foam is hook & loop sensitive allowing maximum adjustability
- Two buckles on rear felt pad are used for tension adjustment

18



BACKS



TWIN V® UNIVERSAL


- Two Universal sizes (S-L, XL-XXXL ) reduce inventory
- Lordotic Molding provides support and intimate fit
- Intimate fit prevents brace migration, increases patient comfort and compliance
- Patient's control and support is easily adjusted as condition changes



| ORDER No. | DESCRIPTION |
|-----------|-------------|
|-----------|-------------|

|     |  |
|-----|--|
| 482 | TWIN V® UNIVERSAL W/ SAGITTAL PANELS<br>PDAC L0637/0650 APPROVED |
| 481 | TWIN V® UNIVERSAL<br>PDAC L0631/L0648 APPROVED                   |

| SIZING: |           |                 |
|---------|-----------|-----------------|
| S/L     | 30" - 43" | (76CM - 109CM)  |
| XL/XXXL | 43" - 58" | (109CM - 147CM) |

VIEW APPLICATION VIDEO AT HELY-WEBER.COM



TWIN V® UNIVERSAL LITE

- Thermoplastic panels provide optimum anterior and posterior support
- Designed for easy patient application and adjustment with patient assist straps and thumb loops
- Flexible vertical stays provide additional support and step-down options
- Cotton spandex and spacer fabric construction provides cool comfortable support
- Brace may be practitioner modified at the unique trim-line for size adjustments and step-down
- Tapered 8 inch belt to maximize comfort



| ORDER No. | DESCRIPTION |
|-----------|-------------|
|-----------|-------------|

|     |   |
|-----|---|
| 480 | TWIN V® UNIVERSAL LITE<br>PDAC L0627/L0642 APPROVED |
|-----|---|

| SIZING:                          |           |                 |
|----------------------------------|-----------|-----------------|
| MEASURE WIDEST PART OF HIP/WAIST |           |                 |
| S/L                              | 30" - 43" | (76CM - 109CM)  |
| XL/XXXL                          | 43" - 58" | (109CM - 147CM) |



ORDER No. DESCRIPTION

- 418

PDAC

L0650

APPROVED

TWIN V® W/ SAGITTAL PANELS
- 408

PDAC

L0648

APPROVED

TWIN V®

| SIZING: | MEASURE WIDEST PART OF HIP/WAIST |
|---------|----------------------------------|
| S       | 27" - 32" (69CM - 81CM)          |
| M       | 33" - 38" (84CM - 97CM)          |
| L       | 37" - 42" (94CM - 107CM)         |
| XL      | 41" - 46" (104CM - 117CM)        |
| XXL     | 46" - 52" (117CM - 132CM)        |
| XXXL    | 51" - 58" (130CM - 147CM)        |
| 4XL     | 58" - 65" (147CM - 165CM)        |



VIEW APPLICATION VIDEO AT HELY-WEBER.COM



TWIN V® SPINAL ORTHOSIS

- Thermoplastic panels provide optimum anterior and posterior support
- Two tension straps provide adjustable intracavitary pressure
- Flexible vertical spiral stays provide focused support
- Patient assist straps and finger assist loops for easy application
- Breathable cotton spandex and spacer fabric construction keeps the patient cool and comfortable
- Available with sagittal panels for additional lateral support

ORDER No. DESCRIPTION

- 419

PDAC

L0650

APPROVED

TWIN V® LOW PROFILE W/ SAGITTAL PANELS
- 409

PDAC

L0648

APPROVED

TWIN V® LOW PROFILE

| SIZING: | MEASURE WIDEST PART OF HIP/WAIST |
|---------|----------------------------------|
| S       | 27" - 32" (69CM - 81CM)          |
| M       | 33" - 38" (84CM - 97CM)          |
| L       | 37" - 42" (94CM - 107CM)         |
| XL      | 41" - 46" (104CM - 117CM)        |
| XXL     | 46" - 52" (117CM - 132CM)        |
| XXXL    | 51" - 58" (130CM - 147CM)        |



VIEW APPLICATION VIDEO AT HELY-WEBER.COM



TWIN V® LOW PROFILE

- Thermoplastic panels provide optimum anterior and posterior support
- Two tension straps provide adjustable intracavitary pressure
- Flexible vertical spiral stays provide focused support
- Patient assist straps and finger assist loops for easy application
- Breathable cotton spandex and spacer fabric construction keeps the patient cool and comfortable
- Available with sagittal panels for additional lateral support
- Tapered 8 inch belt to maximize comfort

ORDER No. DESCRIPTION

- 420

PDAC

L0642

APPROVED

TWIN V® LITE

| SIZING: | MEASURE WIDEST PART OF HIP/WAIST |
|---------|----------------------------------|
| S       | 27" - 32" (69CM - 81CM)          |
| M       | 33" - 38" (84CM - 97CM)          |
| L       | 37" - 42" (94CM - 107CM)         |
| XL      | 41" - 46" (104CM - 117CM)        |
| XXL     | 46" - 52" (117CM - 132CM)        |
| XXXL    | 51" - 58" (130CM - 147CM)        |



TWIN V® LITE

- Thermoplastic panels provide optimum anterior and posterior support
- Designed for easy patient application and adjustment with patient assist straps and thumb loops
- Flexible vertical spiral stays provide focused support
- Cotton spandex and spacer fabric construction provides cool comfortable support
- Brace may be practitioner modified at the unique trim-line for size adjustments
- Tapered 8 inch belt to maximize comfort

LEI BACK™ II

- Unique design available in Kuhl™ Neoprene
- A posterior pocket is provided for the lordotic foam pad for step-down support
- Anterior and posterior panel contours for a personalized anatomical fit



ORDER No. DESCRIPTION

- 5805

PDAC

L0642

APPROVED

KUHL™ LEI BACK™ II

| SIZING: | MEASURE WIDEST PART OF HIP/WAIST |
|---------|----------------------------------|
| XS      | 24" - 34" (61CM - 86CM)          |
| REG     | 34" - 46" (86CM - 117CM)         |
| XL      | 45" - 58" (114CM - 147CM)        |

LEI BACK™ LUMBOSACRAL BRACE

- Kuhl™ perforated neoprene construction provides comfortable compression
- 4-way stretch of neoprene conforms to a variety of body contours
- Double pull provides support and compression
- 4 flexible spiral stays prevent brace rollover



ORDER No. DESCRIPTION

- 5809

PDAC

L0625

APPROVED

KUHL™ LEI BACK™ BRACE

| SIZING: | MEASURE WIDEST PART OF HIP/WAIST |
|---------|----------------------------------|
| XS      | 24" - 34" (61CM - 86CM)          |
| REG     | 34" - 46" (86CM - 117CM)         |
| XL      | 45" - 58" (114CM - 147CM)        |

DOUBLE BACK LSO™

- Double side pulls with patient assist straps for easy application
- Heavy elastic construction ensures comfort
- Elastic pocket accommodates polypro insert for support
- Supports the lower back and encourages normal range of motion



ORDER No. DESCRIPTION

- 103

PDAC

L0641

APPROVED

DOUBLE BACK LSO™ W/ POSTERIOR PLASTIC INSERT

| SIZING: | MEASURE WIDEST PART OF HIP/WAIST |
|---------|----------------------------------|
| XS      | 26" - 30" (66CM - 76CM)          |
| S       | 30" - 34" (76CM - 86CM)          |
| M       | 34" - 38" (86CM - 97CM)          |
| L       | 38" - 42" (97CM - 107CM)         |
| XL      | 42" - 46" (107CM - 117CM)        |
| XXL     | 46" - 50" (117CM - 127CM)        |
| XXXL    | 50" - 54" (127CM - 137CM)        |



VIEW APPLICATION VIDEO AT HELY-WEBER.COM



LOWER  
EXTREMITY



GLOBAL KNEE® AC

- Adjustable Calf allows brace to fit even the hard to fit “conical legged” patients
- Provides stability and addresses patellar pain associated with early and end stage OA patients
- Complex construction of spacer fabrics, Lycra, and urethane laminate creates a lightweight breathable base
- Rigid range of motion hinges block peak varus/valgus movement during ambulation
- Patellar sling with Trakaderm® buttress lifts and elevates the patella into a more natural position for increased stability and less pain
- Provides compression for chronic swelling
- Soft, comfortable three-layer condyle pad system
- The Patella “Sling” allows for directional control of the patella



| ORDER No.                        | DESCRIPTION                |
|----------------------------------|----------------------------|
| 7640-AC                          | GLOBAL KNEE® AC            |
|                                  | PDAC L1832/L1833 APPROVED  |
|                                  | L2795 ADD ON CODE OPTIONAL |
|                                  | USPN 10,206,804            |
| SIZING: MEASURE 6" ABOVE PATELLA |                            |
| XS                               | 15" - 17" (38CM-43CM)      |
| S                                | 16" - 18" (41CM-46CM)      |
| M                                | 17.5" - 19.5" (44CM-50CM)  |
| L                                | 18.5" - 20.5" (47CM-52CM)  |
| XL                               | 19.5" - 22" (50CM-56CM)    |
| XXL                              | 20.5" - 23" (52CM-58CM)    |
| XXXL                             | 24" - 27" (61CM-69CM)      |

LARGER AND CUSTOM SIZES FOR 7640  
MODELS AVAILABLE UPON REQUEST



VIEW APPLICATION VIDEO  
AT HELY-WEBER.COM





ORDER No. DESCRIPTION GLOBAL KNEE®

7640 GLOBAL KNEE®  
PDAC L1832/L1833 APPROVED  
L2795 ADD ON CODE OPTIONAL  
USPN 10,206,804

| SIZING: USE LARGEST OF TWO MEAS. |   |                                      |
|----------------------------------|---|--------------------------------------|
|                                  | CALF<br>MEASURE AT WIDEST<br>PART OF CALF | THIGH<br>MEASURE 6" ABOVE<br>PATELLA |
| XS                               | 11"-13"<br>(28CM-33CM)                    | 15"-17"<br>(38CM-43CM)               |
| S                                | 12.5"-13.5"<br>(32CM-34CM)                | 16"-18"<br>(41CM-46CM)               |
| M                                | 13.5"-15.5"<br>(34CM-39CM)                | 17.5"-19.5"<br>(44CM-50CM)           |
| L                                | 14.5"-16.5"<br>(37CM-42CM)                | 18.5"-20.5"<br>(47CM-52CM)           |
| XL                               | 16.5"-19"<br>(42CM-48CM)                  | 19.5"-22"<br>(50CM-56CM)             |
| XXL                              | 18"-20.5"<br>(46CM-52CM)                  | 20.5"-23"<br>(52CM-58CM)             |
| XXXL                             | 19"-22"<br>(48CM-56CM)                    | 24"-27"<br>(61CM-69CM)               |

LARGER AND CUSTOM SIZES FOR  
7640 MODELS AVAILABLE UPON REQUEST



VIEW APPLICATION VIDEO  
AT HELY-WEBER.COM



- Hybrid design of a fixed calf cuff and anterior closure thigh makes it easy to apply
- Provides stability and addresses patellar pain associated with early and end stage OA patients
- Complex construction of spacer fabrics, Lycra and urethane laminate creates a lightweight breathable base
- Rigid range of motion hinges block peak varus/valgus movement during ambulation
- Patellar sling with Trakaderm® buttress lifts and elevates the patella into a more natural position for increased stability and less pain
- Provides compression for chronic swelling
- Soft, comfortable three layer condyle pad system

ORDER No. DESCRIPTION GK ORIGIN

5640 GK ORIGIN  
PDAC L1820 APPROVED  
USPN 10,206,804

| SIZING: MEASURE CIRCUMFERENCE 6" ABOVE MID-PATELLA WITH LEG FULLY EXTENDED |           |               |
|--|-----------|---------------|
| XS   | 13" - 15" | (33CM - 38CM) |
| S  | 15" - 17" | (38CM - 43CM) |
| M  | 17" - 20" | (43CM - 51CM) |
| L  | 19" - 22" | (48CM - 56CM) |
| XL   | 21" - 24" | (53CM - 61CM) |
| XXL  | 23" - 27" | (58CM - 69CM) |
| XXXL   | 26" - 29" | (66CM - 74CM) |



THE ADJUSTABLE  
CALF ALLOWS THE  
BRACE TO BE FIT  
ON A CONICAL  
SHAPED LEG

- Adjustable calf allows brace to fit even the hard to fit patients
- Provides stability and addresses patellar pain associated with early and end stage OA patients
- Rigid polycentric hinges block peak varus/valgus movement during ambulation
- Patellar sling lifts and elevates the patella into a more natural position for increased stability and less pain
- Provides compression for chronic swelling

KNAPP™II ORDER No. DESCRIPTION

- Anterior closure design allows the brace to be contoured to give a custom fit
- Combinations of spacer mesh and double knit material provides compression support and edema control
- Low profile hinge
- Sleek design keeps hinge close to the knee for optimum support
- Pop Stop™ ROM system (patent pending) allows setting changes in seconds
- Additional ROM stops included with brace
- Condyle pads for added comfort
- Lightweight and breathable material forms to the leg wicks moisture away from body as well as keeping the leg cool during activities
- Malleable uprights allow the uprights to be bent to the contour of the patient for maximum support and comfort
- 12" knee brace ideal for MCL/LCL sprains and strains
- 16" knee brace ideal for ACL/PCL, MCL/LCL sprains and strains



KNAPP™II  
HINGED KNEE ORTHOSIS  
PDAC L1832/L1833 APPROVED  
7686 12" KNAPP™II HINGED KNEE  
7688 16" KNAPP™II HINGED KNEE  
PATENT PENDING

| SIZING: MEASUREMENTS ARE TAKEN AT 6" ABOVE MID-PATELLA |           |            |
|--|-----------|------------|
| XS   | 13" - 15" | 33 - 38 CM |
| S  | 15" - 18" | 38 - 46 CM |
| M  | 18" - 20" | 46 - 50 CM |
| L  | 20" - 23" | 50 - 58 CM |
| XL   | 23" - 25" | 58 - 64 CM |
| XXL  | 25" - 28" | 64 - 71 CM |

LARGER AND CUSTOM SIZES FOR  
7686/7688 MODELS AVAILABLE  
UPON REQUEST

GREATER ROM LIMITS AVAILABLE IN  
PROFESSIONAL STOP KIT #X1224-K



VIEW APPLICATION VIDEO  
AT HELY-WEBER.COM

THOMAS P. KNAPP, M.D.  
Cedars-Sinai®  
Kerlan-Jobe Institute



7686

7688



ORDER No. DESCRIPTION KNAPP™ HINGED KNEE ORTHOSIS - ANTERIOR CLOSURE

KNAPP™ HINGED KNEE ORTHOSIS *KUHL™* NEOPREN  
5656 12" STANDARD HINGES  
5658 16" STANDARD HINGES  
5656HH 12" HEAVY HINGES  
5658HH 16" HEAVY HINGES

PDAC L1820 APPROVED

5656ROM 12" RANGE OF MOTION  
5658ROM 16" RANGE OF MOTION

PDAC L1832/L1833 APPROVED

NOSWET® KNAPP™ HINGED KNEE ORTHOSIS  
7656 12" STANDARD HINGES  
7658 16" STANDARD HINGES

7656HH 12" HEAVY HINGES  
7658HH 16" HEAVY HINGES

PDAC L1820 APPROVED

7656ROM 12" RANGE OF MOTION

PDAC L1832/L1833 APPROVED

SIZING: MEASURE AT MID-PATELLA WITH LEG FULLY EXTENDED

|     |           |               |
|-----|-----------|---------------|
| XS  | 10" - 12" | (25CM - 31CM) |
| S   | 12" - 14" | (31CM - 36CM) |
| M   | 14" - 16" | (36CM - 41CM) |
| L   | 16" - 18" | (41CM - 46CM) |
| XL  | 18" - 20" | (46CM - 51CM) |
| XXL | 20" - 22" | (51CM - 56CM) |

LARGER & CUSTOM SIZES FOR 5656/5658  
MODELS AVAILABLE UPON REQUEST



THOMAS P. KNAPP, M.D.  
Cedars-Sinai®  
Kerlan-Jobe Institute



5656/7656

- Available in *Kuhl™* neoprene or *NoSwet®* non-neoprene material
- 12" & 16" length (30.5cm & 40.6cm)
- 3 hinge styles
- Easily accessible hinge allows interchangeability of hinge styles
- Hinges are fully covered for sports activity
- Elasticized thigh strap ensures comfort
- Ideal for pre/post-operative care of MCL/LCL and meniscal repairs
- Longer length can be used for mild to moderate ACL/PCL sprains & strains and as an adjunct to patella fracture



5658/7658

ORDER No. DESCRIPTION VELOCITY HINGED KNEE - PULL-UP/ANTERIOR CLOSURE

5645 VELOCITY HINGED KNEE  
5647 12" PULL-UP  
5646 16" PULL-UP  
5648 12" ANTERIOR CLOSURE  
5648 16" ANTERIOR CLOSURE

PDAC L1832/L1833 APPROVED

SIZING MEASURE CIRCUMFERENCE 6"  
ABOVE MID-PATELLA WITH  
LEG FULLY EXTENDED

| VELOCITY HINGED KNEE PULL-UP |           |               |
|------------------------------|-----------|---------------|
| XS                           | 13" - 15" | (33CM - 38CM) |
| S                            | 15" - 17" | (38CM - 43CM) |
| M                            | 16" - 18" | (41CM - 46CM) |
| L                            | 18" - 20" | (46CM - 51CM) |
| XL                           | 20" - 22" | (51CM - 56CM) |
| XXL                          | 22" - 24" | (56CM - 61CM) |
| XXXL                         | 24" - 26" | (61CM - 66CM) |

VELOCITY HINGED KNEE ANTERIOR CLOSURE

|      |           |               |
|------|-----------|---------------|
| XS   | 12" - 16" | (31CM - 41CM) |
| S    | 14" - 18" | (36CM - 46CM) |
| M    | 16" - 20" | (41CM - 51CM) |
| L    | 17" - 22" | (43CM - 56CM) |
| XL   | 19" - 25" | (48CM - 64CM) |
| XXL  | 23" - 28" | (58CM - 71CM) |
| XXXL | 25" - 29" | (63CM - 74CM) |



5646

- Constructed of *Kuhl™* perforated neoprene to maximize patient comfort and compliance
- Range-of-Motion hinges and stop kit allows for injury specific ROM settings
- Brace can be locked at full extension
- Indicated and ideally suited for treatment of ACL, PCL, MCL and LCL sprains, as well as post-operative and non-operative meniscus injuries
- Available in anterior closure and pull-up designs
- Removable hinges for easy brace cleaning
- Interchangeable ext/flex stop provide 10 degrees increases in control



5648



5645



5647

SHIELDS II® BRACE

- Provides dynamic control of the patella reducing discomfort & enhancing recovery
- Trakaderm® U-shaped buttress allows for treatment of multiple patella indications and fine tuning the position of the patella no matter the malalignment
- One brace fits both left & right for less Inventory
- Low Profile and extremely durable
- Independent calf and thigh straps adjust for contours of the leg
- Hinges provide medial and lateral stability can be removed for step-down application (6675)



6674



6675

ORDER No. DESCRIPTION

6674 SHIELDS II® BRACE  
6675 HINGED SHIELDS II® BRACE  
PDAC L1820 APPROVED

| SIZING: | MEASURE 4" ABOVE MID-PATELLA |               |
|---------|------------------------------|---------------|
| XS      | 10" - 11.5"                  | (25CM - 29CM) |
| S       | 12" - 13.5"                  | (30CM - 34CM) |
| M       | 14" - 16"                    | (36CM - 41CM) |
| L       | 16.5" - 19"                  | (42CM - 48CM) |
| XL      | 19.5" - 22"                  | (50CM - 56CM) |
| XXL     | 22.5" - 25"                  | (57CM - 64CM) |

| MEASURE AT MID-PATELLA |             |               |
|------------------------|-------------|---------------|
| XS                     | 10" - 11"   | (25CM - 28CM) |
| S                      | 11.5" - 13" | (29CM - 33CM) |
| M                      | 13.5" - 15" | (34CM - 38CM) |
| L                      | 15" - 17"   | (38CM - 43CM) |
| XL                     | 17" - 19"   | (43CM - 48CM) |
| XXL                    | 19" - 21"   | (48CM - 53CM) |



VIEW APPLICATION VIDEO  
AT HELY-WEBER.COM

CLARENCE L. SHIELDS, M.D.  
Cedars-Sinai®  
Kerlan-Jobe Institute



KNEESIO™ BRACE

- Extremely light weight "only 8 oz."
- Flexible hinge contours to patients leg
- Breathable material keeps athlete cool while stabilizing patella tracking
- Four adjustable patella straps provides customized patellar support
- Removable/adjustable buttress for additional control where needed



ORDER No. DESCRIPTION

7675 KNEESIO™ BRACE  
PDAC L1820 APPROVED

| SIZING: | MEASURE 6" ABOVE MID-PATELLA |               |
|---------|------------------------------|---------------|
| XS      | 13" - 15"                    | (33CM - 38CM) |
| S       | 15" - 18"                    | (38CM - 46CM) |
| M       | 18" - 20"                    | (46CM - 50CM) |
| L       | 20" - 23"                    | (50CM - 58CM) |
| XL      | 23" - 25"                    | (58CM - 64CM) |
| XXL     | 25" - 28"                    | (64CM - 71CM) |



VIEW APPLICATION VIDEO  
AT HELY-WEBER.COM

WENDELL BECTON, M.D.  
The Orthopedic Center  
of St. Louis





ORDER No.

DESCRIPTION

5675

SHIELDS® HINGED PATELLA STABILIZER (ORIGINAL)

5674

KUHL™ SHIELDS® BRACE (ORIGINAL)

SIZING:

MEASURE CIRCUMFERENCE 6" ABOVE MID-PATELLA WITH LEG FULLY EXTENDED

|     |           |               |
|-----|-----------|---------------|
| XS  | 10" - 12" | (25CM - 30CM) |
| S   | 12" - 14" | (30CM - 36CM) |
| M   | 14" - 16" | (36CM - 41CM) |
| L   | 16" - 18" | (41CM - 46CM) |
| XL  | 18" - 20" | (46CM - 51CM) |
| XXL | 20" - 24" | (51CM - 61CM) |



VIEW APPLICATION VIDEO AT HELY-WEBER.COM

CLARENCE L. SHIELDS, M.D.  
Cedars-Sinai®  
Kerlan-Jobe Institute

SHIELDS® BRACE (ORIGINAL)

- Provides dynamic control of the patella reducing discomfort and enhancing recovery
- Trakaderm® U-shaped buttress allows for treatment of multiple patella indications
- The directional force strap and the Trakaderm® buttress allows for fine tuning of patella support and maintains the patella in the desired position
- One brace fits both left & right for less inventory
- Independent calf and thigh straps adjust for contours of the leg
- Hinges provide medial and lateral stability and can be removed for step-down application (5675)



ORDER No.

DESCRIPTION

5645PS

12" VELOCITY HINGED KNEE RANGE OF MOTION W/ PATELLA STABILIZATION - PULL UP

5646PS

12" VELOCITY HINGED KNEE RANGE OF MOTION W/ PATELLA STABILIZATION - ANTERIOR CLOSURE

PDAC L1832, L1833 APPROVED

L2795 ADD ON CODE OPTIONAL

SIZING:

MEASURE CIRCUMFERENCE 6" ABOVE MID-PATELLA WITH LEG FULLY EXTENDED

|      | PULL UP   | ANTERIOR CLOSURE |
|------|-----------|------------------|
| XS   | 13" - 15" | 12" - 16"        |
| S    | 15" - 17" | 14" - 18"        |
| M    | 16" - 18" | 16" - 20"        |
| L    | 18" - 20" | 17" - 22"        |
| XL   | 20" - 22" | 19" - 25"        |
| XXL  | 22" - 24" | 23" - 28"        |
| XXXL | 24" - 26" | 25" - 29"        |



VIEW APPLICATION VIDEO AT HELY-WEBER.COM

VELOCITY PS PULL UP/ANTERIOR CLOSURE

- The brace offers MCL and LCL support and patella control
- Universal Left/Right U-shaped Trakaderm® buttress can be moved to opposite side by simply disconnecting the hook and loop flap and applying the device to the opposite upright
- Kuhl™ perforated neoprene breathes and offers warmth & compression
- Range of Motion hinges feature a polycentric adjustable design to control extension and flexion
- Circumferential straps have elastic hamstring tension relief
- Available in pull-up and anterior closure
- Interchangeable ext/flex stop provides 10 degrees increase in control



ORDER No.

DESCRIPTION

3655

12" (30.5CM) AXIS HINGED KNEE SLEEVE - OPEN POPLITEAL

PDAC L1812 APPROVED

SIZING:

MEASURE AT MID-PATELLA WITH LEG FULLY EXTENDED

|     |           |               |
|-----|-----------|---------------|
| XS  | 12" - 13" | (30CM - 33CM) |
| S   | 13" - 14" | (33CM - 36CM) |
| M   | 14" - 15" | (36CM - 38CM) |
| L   | 15" - 17" | (38CM - 43CM) |
| XL  | 16" - 18" | (41CM - 46CM) |
| XXL | 18" - 20" | (46CM - 51CM) |

LARGER AND CUSTOM SIZES ARE AVAILABLE UPON REQUEST

AXIS HINGED KNEE SLEEVE

- 1/8" (3mm) Neoprene with nylon outer covering
- Polycentric hinges prevent hyperextension
- Covered hinges allow full sports participation
- Pull up loops aid in application
- Circumferential straps have elastic hamstring tension relief



ORDER No.

DESCRIPTION

24100

BASIC KNEE SLEEVE CLOSED PATELLA, CLOSED POPLITEAL

3605

OPEN PATELLA, OPEN POPLITEAL

24105

OPEN PATELLA, CLOSED POPLITEAL

5607

KUHL™ KNAPP™ SLEEVE ANTERIOR CLOSURE

SIZING:

MEASURE AT MID-PATELLA W/ LEG FULLY EXTENDED

|     |           |               |
|-----|-----------|---------------|
| XS  | 12" - 13" | (30CM - 33CM) |
| S   | 13" - 14" | (33CM - 36CM) |
| M   | 14" - 15" | (36CM - 38CM) |
| L   | 15" - 17" | (38CM - 43CM) |
| XL  | 16" - 18" | (41CM - 46CM) |
| XXL | 18" - 20" | (46CM - 51CM) |

LARGER AND CUSTOM SIZES AVAILABLE UPON REQUEST

BASIC KNEE SLEEVE/KNAPP™ SLEEVE ANTERIOR CLOSURE

- Basic Knee Sleeve**
  - 1/8" (3mm) Neoprene with 2-sided nylon outer covering
  - Provides compression and warmth
- Knapp™ Sleeve - Anterior Closure**
  - 1/8" (3mm) Perforated neoprene construction for cool comfortable compression
  - Easy application, even over clothes
  - Allows extra adjustment for sizing and use with dressings



ORDER No.

DESCRIPTION

203

KNEE IMMOBILIZER ELASTIC STRAPS 16" (41CM)

204

18" (46CM)

205

20" (51CM)

206

22" (56CM)

SIZING:

FITS MID-PATELLA CIRC. 12.0"-23.0" (30CM - 58CM)

KNEE IMMOBILIZER - ELASTIC STRAPS

- Made from Trimable™ 3/8" (9.5mm) pile foam laminate
- Reinforced, moldable posterior stays provide additional rigidity
- Medial and lateral movable stay pods offer extra support
- Contact closure application straps
- Elastic adjustment straps provide added compression
- Pile exterior is hook & loop sensitive for easy adjustment



ORDER No.

DESCRIPTION / SIZE

211

3 PANEL KNEE IMMOBILIZER 12" (30CM)

213

16" (41CM)

214

18" (46CM)

215

20" (51CM)

216

22" (56CM)

217

24" (61CM)

SIZING:

FITS MID-PATELLA CIRC. 12" - 24" (30CM - 61CM)

THREE PANEL KNEE IMMOBILIZER

- Provides rigid support for the knee
- 3/8" (9.5mm) pile foam laminate with black flannel interior ensures patient comfort
- Two detachable panels for circumferential adjustments for custom fitting
- Posterior and medial/lateral stays are malleable
- Optional patella strap provides extra knee control
- Available in 6 lengths





ORDER No.

DESCRIPTION

HINGED LATERAL J PATELLA STABILIZER

5694H

KUHL™ HINGED LATERAL J STABILIZER - OPEN POPLITEAL WITH CONDYLE PADS

PDAC

L1820

APPROVED

SIZING:

MEASURE AT MID-PATELLA WITH LEG FULLY EXTENDED

|     |           |               |
|-----|-----------|---------------|
| XS  | 12" - 13" | (30CM - 33CM) |
| S   | 13" - 14" | (33CM - 36CM) |
| M   | 14" - 15" | (36CM - 38CM) |
| L   | 15" - 17" | (38CM - 43CM) |
| XL  | 16" - 18" | (41CM - 46CM) |
| XXL | 18" - 20" | (46CM - 51CM) |

LARGER AND CUSTOM SIZES ARE AVAILABLE UPON REQUEST

SPECIFY RIGHT OR LEFT



- J Buttress in conjunction with Anterior Flap gives dynamic patella control
- Biaxial Hinges provide medial and lateral stability
- Kuhl™ perforated neoprene allows moisture to escape and the skin to remain dry

ORDER No.

DESCRIPTION

UNIVERSAL BUTTRESS HINGED PATELLA STABILIZER

3630H

UNIVERSAL BUTTRESS HINGED PATELLA STABILIZER, OPEN POPLITEAL

PDAC

L1812

APPROVED

SIZING:

MEASURE AT MID-PATELLA WITH LEG FULLY EXTENDED

|     |           |               |
|-----|-----------|---------------|
| XS  | 12" - 13" | (30CM - 33CM) |
| S   | 13" - 14" | (33CM - 36CM) |
| M   | 14" - 15" | (36CM - 38CM) |
| L   | 15" - 17" | (38CM - 43CM) |
| XL  | 16" - 18" | (41CM - 46CM) |
| XXL | 18" - 20" | (46CM - 51CM) |

LARGER AND CUSTOM SIZES ARE AVAILABLE UPON REQUEST



- Universal tubular buttress is movable/removable for individualized patient care
- 12" (30.5cm) Nylon covered neoprene base offers comfort and compression
- Wide bi-axial hinges provide medial/lateral ligamentous support
- Pull up loops aid in application
- Neoprene hinge covers allow for full sports participation

ORDER No.

DESCRIPTION

INFERIOR U BUTTRESS HINGED PATELLA STABILIZER

3641H

INFERIOR U BUTTRESS HINGED PATELLA STABILIZER, OPEN POPLITEAL

PDAC

L1812

APPROVED

SIZING:

MEASURE AT MID-PATELLA WITH LEG FULLY EXTENDED

|     |           |               |
|-----|-----------|---------------|
| XS  | 12" - 13" | (30CM - 33CM) |
| S   | 13" - 14" | (33CM - 36CM) |
| M   | 14" - 15" | (36CM - 38CM) |
| L   | 15" - 17" | (38CM - 43CM) |
| XL  | 16" - 18" | (41CM - 46CM) |
| XXL | 18" - 20" | (46CM - 51CM) |

LARGER AND CUSTOM SIZES ARE AVAILABLE UPON REQUEST



- Inferior tubular buttress provides patella control
- 12" (30.5cm) Nylon covered neoprene offers comfort and compression
- Wide bi-axial hinges provide medial/lateral ligamentous support
- Neoprene hinge covers allow full sports participation

ORDER No.

DESCRIPTION

MATT STRAP™ II UNIVERSAL

6709

MATT STRAP™ II UNIVERSAL


SIZING:


ENCIRCLE THE AREA OF APPLICATION TO MEASURE

|            |             |               |
|------------|-------------|---------------|
| FITS CIRC. | OF 7" - 20" | (18CM - 51CM) |
|------------|-------------|---------------|

MATTHEW SHAPIRO, M.D.

SLOCUM CENTER FOR ORTHOPEDICS & SPORTS MEDICINE





- Malleable aluminum bolster conforms to the contours of the tibia for improved function and comfort
- The strap and buttress combination reduces adverse stress on the patella tendon
- Shock absorbing buttress diffuses energy and reduces the load on the patella tendon
- Pain reduction allows the patient to resume activity

ORDER No.

DESCRIPTION

MATT STRAP™ - SHOCK ABSORPTION

3704

MATT STRAP™

SIZING:

MEASURE CIRC. OF AFFECTED AREA

|     |             |               |
|-----|-------------|---------------|
| XS  | 7" - 8"     | (18CM - 20CM) |
| S   | 8" - 9"     | (20CM - 23CM) |
| M   | 9" - 10.5"  | (23CM - 27CM) |
| L   | 10.5" - 12" | (27CM - 30CM) |
| XL  | 12" - 14"   | (30CM - 36CM) |
| XXL | 13" - 16"   | (33CM - 41CM) |



- Useful for patella tendonitis, "Jumpers Knee", and minor knee pain
- Loop & lock design provides focused compression over the shock absorbing buttress

ORDER No.

DESCRIPTION

SEMI-UNIVERSAL MATT STRAP™

3707 (S/M)  
3707 (L/XL)

SEMI-UNIVERSAL MATT STRAP™  
SEMI-UNIVERSAL MATT STRAP™

SIZING:

MEASURE CIRC. OF AFFECTED AREA

|      |           |               |
|------|-----------|---------------|
| S/M  | 11" - 14" | (28CM - 36CM) |
| L/XL | 14" - 16" | (36CM - 41CM) |



- Useful for patella tendonitis, "Jumpers Knee", and minor knee pain
- Non-loop lock construction
- Unique double pull compression straps
- Two sizes
- The elastic double-pull strap provides unique additional compression and adjustment

ORDER No.

DESCRIPTION

WEBLY® ANKLE ORTHOSIS

304

PDAC

L1902

APPROVED

WEBLY® ANKLE ORTHOSIS

308A

OPTIONAL STAYS (PAIR)

SIZING:

FITTING BY SHOE SIZE

|       |        |           |      |
|-------|--------|-----------|------|
| MEN   | NARROW | REGULAR   | WIDE |
| XXS   | 6.5    | 5.5 - 6   | 5    |
| XS    | 7.5    | 6.5 - 7   | 6    |
| S     | 8.5    | 7.5 - 8   | 7    |
| M     | 11.5   | 8.5 - 11  | 8    |
| L     | 13.5   | 11.5 - 13 | 10   |
| XL    | 14.5   | 13.5 - 14 | 12   |
| XXL   |        | 14.5 - 16 | 13   |
| WOMEN | NARROW | REGULAR   | WIDE |
| XXS   | 7.5    | 6.5 - 7   | 6    |
| XS    | 9      | 7.5 - 8.5 | 7    |
| S     | 10.5   | 8 - 10    | 7.5  |
| M     | 12.5   | 10.5 - 12 | 10   |
|       | 14.5   | 12 - 14   | 11   |



- Constructed of heavy gauge ballistic nylon, with non-stretch, non-fraying nylon straps
- Limits inversion/eversion while not restricting plantar and dorsiflexion
- Neoprene tongue and Achilles relief padding aid in player/patient comfort
- Fits in shoes and work boots
- Seamless design maximizes comfort

ORDER No.

DESCRIPTION

318

PDAC

L1902

APPROVED

RAPID ZAP™ ANKLE ORTHOSIS

308A

OPTIONAL STAYS (PAIR)

SIZING:

FITTING BY SHOE SIZE

| MEN | NARROW | REGULAR   | WIDE |
|-----|--------|-----------|------|
| XXS | 6.5    | 5.5 - 6   | 5    |
| XS  | 7.5    | 6.5 - 7   | 6    |
| S   | 8.5    | 7.5 - 8   | 7    |
| M   | 11.5   | 8.5 - 11  | 8    |
| L   | 13.5   | 11.5 - 13 | 10   |
| XL  | 14.5   | 13.5 - 14 | 12   |
| XXL |        | 14.5 - 16 | 13   |

| WOMEN | NARROW | REGULAR   | WIDE |
|-------|--------|-----------|------|
| XXS   | 7.5    | 6.5 - 7   | 6    |
| XS    | 9      | 7.5 - 8.5 | 7    |
| S     | 10.5   | 8 - 10    | 7.5  |
| M     | 12.5   | 10.5 - 12 | 10   |
| L     | 14.5   | 12 - 14   | 11   |



RAPID ZAP™ ANKLE ORTHOSIS

- Unique lacing design allows independent forefoot, midfoot and ankle adjustments
- Figure 8 straps mimic ankle taping in a preferred neutral position providing resistance to “roll-over” ankle injury while not restricting walking and running
- Sturdy construction and seamless design offers comfort, application, quality and durability
- Elastic cuff closure secures laces and straps while offering additional compression
- Adjustable lacing panel enables the opening to accommodate wide or swollen ankles
- Seamless design maximizes comfort
- Quick on and off application

ORDER No.

DESCRIPTION

318H

PDAC

L1906

APPROVED

HINGED RAPID ZAP™

SIZING:

FITTING BY SHOE SIZE

| WOMEN     | NARROW | REGULAR | WIDE |
|-----------|--------|---------|------|
| 5 - 6     |        | XS      | S    |
| 6.5 - 8   | XS     | S       | M    |
| 8.5 - 10  | S      | M       | L    |
| 10.5 - 12 | M      | M       | L    |
| 12.5 - 14 | L      | L       | XL   |

| MEN       | NARROW | REGULAR | WIDE |
|-----------|--------|---------|------|
| 6 - 7     |        | XS      | S    |
| 7.5 - 9   | XS     | S       | M    |
| 9.5 - 11  | S      | M       | L    |
| 11.5 - 13 | M      | L       | XL   |
| 13.5 - 15 | L      | XL      | XXL  |
| 15 +      | XL     | XXL     | XXL  |



HINGED RAPID ZAP™

- Sturdy ballistic nylon assembly provides durability and longevity
- Seamless construction eliminates potential irritation
- Low profile hinge provides rigidity and support to the ankle
- Close fitting design slips easily into the shoe
- Large throat opens for long feet and large ankles. Get the correct brace from the beginning. No need to upsize for a swollen ankle
- Universal left/right hinge and footplate assembly extends from midfoot to hind-foot for maximum control and comfort
- Hinged assembly accommodates patients' anatomy
- Figure 8 straps mimic ankle taping to secure the ankle mortise and control inversion/eversion
- Three zone speed lacer design, no need to tighten the laces at each eyelet
- Three speed lacing tabs allow fast independent adjustments to the forefoot, midfoot and ankle

VIEW APPLICATION VIDEO AT HELY-WEBER.COM

ORDER No.

DESCRIPTION

314-BLK

PDAC

L1902

APPROVED

WEBLY ZAP® ANKLE ORTHOSIS

308A

OPTIONAL STAYS (PAIR)

SIZING:

FITTING BY SHOE SIZE

| MEN | NARROW | REGULAR   | WIDE |
|-----|--------|-----------|------|
| XXS | 6.5    | 5.5 - 6   | 5    |
| XS  | 7.5    | 6.5 - 7   | 6    |
| S   | 8.5    | 7.5 - 8   | 7    |
| M   | 11.5   | 8.5 - 11  | 8    |
| L   | 13.5   | 11.5 - 13 | 10   |
| XL  | 14.5   | 13.5 - 14 | 12   |
| XXL |        | 14.5 - 16 | 13   |

| WOMEN | NARROW | REGULAR   | WIDE |
|-------|--------|-----------|------|
| XXS   | 7.5    | 6.5 - 7   | 6    |
| XS    | 9      | 7.5 - 8.5 | 7    |
| S     | 10.5   | 8 - 10    | 7.5  |
| M     | 12.5   | 10.5 - 12 | 10   |
| L     | 14.5   | 12 - 14   | 11   |

VIEW APPLICATION VIDEO AT HELY-WEBER.COM

WEBLY ZAP® ANKLE ORTHOSIS

- Limits ankle inversion/eversion while not restricting plantar and dorsiflexion
- Figure 8 stabilizing straps replicate ankle taping
- Low profile brace fits in shoes and work boots
- Useful prophylactically and for ankle sprains
- Seamless design maximizes comfort
- Elastic strap closure improves compression
- Universal fits right or left
- For additional support add stays

ORDER No.

DESCRIPTION

319

PDAC

L1902

APPROVED

WISELOC™

SIZING:

FITTING BY SHOE SIZE

| WOMEN     | NARROW | REGULAR | WIDE |
|-----------|--------|---------|------|
| 5 - 6     | XXS    | XS      | S    |
| 6.5 - 8   | XS     | S       | M    |
| 8.5 - 10  | S      | M       | L    |
| 10.5 - 12 | M      | M       | L    |
| 12.5 - 14 | L      | L       | XL   |

| MEN       | NARROW | REGULAR | WIDE |
|-----------|--------|---------|------|
| 6 - 7     | XXS    | XS      | S    |
| 7.5 - 9   | XS     | S       | M    |
| 9.5 - 11  | S      | M       | L    |
| 11.5 - 13 | M      | L       | XL   |
| 13.5 - 15 | L      | XL      | XXL  |
| 15 +      | XL     | XXL     | XXL  |



WISELOC™

- Three zone speed lacing adjustments provide a fast and customized fit in seconds
- Low profile, easily slips into shoe
- Syndesmosis support system
- Seamless ballistic nylon body design eliminates seams under foot and provides durability and longevity
- Get the correct brace on your patient from the beginning. No need to go up a size for a swollen ankle
- Upper ankle traveler allows the throat to be adjusted for custom fit
- Circumferential ankle pad and adjustable straps provide incremental adjustments for optimum support and comfort

VIEW APPLICATION VIDEO AT HELY-WEBER.COM

BRIAN ENGE  
PT

WALT RUMFELT  
PA-C



ORDER No.

DESCRIPTION

307

AIR ANKLE STIRRUP

PDAC L4350 APPROVED

SIZING:

UNIVERSAL

ADJUST HEEL STRAP FOR COMFORT

FITS UP TO 3" (8CM) WIDTH



- Promotes Edema Reduction and Supportive Protection for acute ankle sprains
- Universal shell design - one size comfortably fits left and right
- Two Air Cell Design aids in active edema reduction
- Heel pad material widths can easily be adjusted
- Thicker heel pad for durability and support
- Air cell pressure is adjustable for better fit

ORDER No.

DESCRIPTION

333

PLANTAR STRETCHING ORTHOSIS

PDAC L4397 APPROVED

SIZING:

SHOE SIZE:

WOMEN'S

US

EUROPE

S

5 - 7

35 - 37.5

M

7 - 10

37.5 - 41

L

10+

41+

MEN'S

US

EUROPE

S

4 - 6

35.5 - 39

M

6 - 9

39 - 42.5

L

9+

42.5+



- Comfortable foam padded foot and calf straps are adjustable with easy close buckles
- Lightweight design and low profile injection molded shell assist in patient comfort and durability
- Soft, removable and washable foam liner
- Slip-resistant tread (not recommended to ambulate in this device)

ORDER No.

DESCRIPTION

327

SILICONE HEEL CUPS

PDAC L3170 APPROVED

SIZING:

FITS MOST CHILDREN'S SHOES

M

FITS MOST WOMEN'S SHOES

L

FITS MOST MEN'S SHOES



- Made from 100% silicone, which helps relieve lower extremity shock and fatigue
- The heel cups are easy to clean, durable and comfortable to wear
- Dual durometer heel spur pad offers pressure relief

ORDER No.

DESCRIPTION

364

MOLDED CAST SHOE

393

WOODEN SOLE POST-OP SHOE

394

FLEXIBLE POST-OP SHOE

SIZING:

SHOE SIZE

364

A

B

PED

7.25" (18CM)

3.7" (9CM)

S

9.4" (24CM)

4.3" (11CM)

M

11" (28CM)

4.9" (12CM)

L

12.5" (32CM)

5.4" (14CM)

393

394

364

364

Diagram of shoe dimensions A and B

SIZING:

SHOE SIZE

393 / 394

WOMENS

MENS

XS

3 - 5

FEMALE SMALL

5 - 7.5

FEMALE MEDIUM

7.5 - 9

7 - 8

FEMALE LARGE

9.5 - 10.5

8.5 - 9.5

MALE SMALL

9.5 - 10.5

8.5 - 9.5

MALE MEDIUM

11 - 12

10 - 11

MALE LARGE

12.5

11.5 - 13.5

X-LARGE

14 - 15

393 NOT AVAILABLE IN XS AND XL

PEDIATRIC

35

36

ORDER No.

DESCRIPTION

450 TINY

TINY TITAN® WRIST

450-XS

TITAN WRIST™ XS

USPN 6,893,410; 6,960,176

L3908 SUGGESTED

SIZING:

TINY

< 4"

(10CM)

XS

4"-6"

(10CM - 15CM)

SPECIFY RIGHT OR LEFT



TINY TITAN® WRIST

- Pediatric splint constructed of soft, comfortable orthopedic felt laminated to brush nylon exterior
- Anatomically contoured palmar stay is malleable for customization and removable for step-down rehabilitation
- Stockinette tongue eases application
- Anatomically designed for a better fit

ORDER No.

DESCRIPTION

455 TINY

TINY TITAN® THUMB

455-XS

TITAN THUMB™ XS

USPN 6,893,410; 6,960,176  
7,033,331; 7,056,298

L3807/3809 SUGGESTED

SIZING:

TINY

< 4"

(10CM)

XS

4"-6"

(10CM - 15CM)

SPECIFY RIGHT OR LEFT



TINY TITAN® THUMB

- Pediatric splint constructed of soft, comfortable orthopedic felt laminated to brush nylon exterior
- Stockinette tongue eases application
- Anatomically designed for a better fit

ORDER No.

DESCRIPTION

5818-PED

KUHL™ MODABBER™ WRIST LONG PEDIATRIC

5819-PED

KUHL™ MODABBER™ WRIST SHORT PEDIATRIC

L3908 SUGGESTED

SIZING:

PEDIATRIC

WRIST CIRCUMFERENCE

4" - 6"

(10CM - 15CM)

SPECIFY RIGHT OR LEFT



MODABBER™ WRIST ORTHOSIS

- Patented patient assist single strap for easy application
- Kuhl™ perforated neoprene construction provides comfortable compression
- 4 way stretch neoprene ensures brace conforms to contour of hand
- Aluminum palmar stay is removable and malleable for a custom fit

TKO® THE KNUCKLE ORTHOSIS

- Provides positioning for fractures or injuries to metacarpals, phalanges, MCP and IP joints
- Four applications in one brace
- Straps are adjustable for fit and comfort (pain & swelling)
- Distal buddy strapping eliminates rotation of the finger
- Padded stays are malleable and allow for customized fit
- Reversible - one splint acts as a Boxer's Fracture splint (ulnar gutter), Radial Gutter splint, and Middle & Ring Finger splint

TINY CAST-AWAY®

- Dual rod stays- provide superior rigidity and support while remaining low profile
- Memory foam olecranon pad - provides comfort and support to prevent any pressure or irritation to vulnerable areas
- Versatile strapping system - allows for customized fit and comfort
- Available in 3 sizes, right/left ambidextrous reducing inventory needs
- Lightweight, low profile design- promotes compliance improving clinical outcomes
- Simple application- promotes good hygiene/skin care and eliminates cast complications and cast saw anxiety



ORDER No.

DESCRIPTION

3848-XS

TKO® (THE KNUCKLE ORTHOSIS)

3848

TKO® (THE KNUCKLE ORTHOSIS)

L3807 SUGGESTED

USPN 7,278,980; 7,442,177;  
7,455,650; 7,402,149;  
7,276,039

SIZING:

X-SMALL

MEASURE CIRCUMFERENCE AT PALMAR CREASE

<7"


(18CM)

REGULAR

7"-10.5"

(18CM-25CM)

SPECIFY RIGHT OR LEFT



ORDER No.

DESCRIPTION

5612

TINY CAST-AWAY®

L3980 / L3982 / L3762 SUGGESTED

SIZING:

SMALL

MEASURE FROM OLECRANON TO ULNAR STYLOID

6.0 - 7.5

(15CM - 19CM)

MEDIUM


7.5 - 8.5

(19CM - 22CM)

LARGE


8.5 - 10.0

(22CM - 25CM)

VIEW APPLICATION VIDEO AT HELY-WEBER.COM

Bryon Izuka, M.D.

Children's Orthopaedics of Hawaii LLC



37

38



ORDER No.

DESCRIPTION

640-S-8-XS

MTC FRACTURE BRACE®

LCODE L3982 SUGGESTED

650-S-8-XS

MTC FRACTURE BRACE® WITH THUMB SPICA

LCODE L3984 SUGGESTED

USPN 9,782,285

SIZING: PLEASE REFERENCE PAGE 14 FOR SIZING



- Controls forearm pronation/supination
- Allows elbow flexion/extension
- Simple and fast application
- Comfortable fit
- Easily adjusted
- Step-down capability

ORDER No.

DESCRIPTION

501-PED

THE GUS SHOULDER IMMOBILIZER SLING

LCODE L3670 SUGGESTED

501-XS

THE GUS SHOULDER IMMOBILIZER SLING

501-S

THE GUS SHOULDER IMMOBILIZER SLING

LCODE L3670 SUGGESTED

SIZING:

MEASURE FROM OLECRANON TO KNUCKLES

HEIGHT

|     |                       |           |
|-----|-----------------------|-----------|
| PED | 7" - 9" (18 - 23CM)   | 5" (12CM) |
| XS  | 8"-10" (20CM - 25CM)  | 7" (18CM) |
| S   | 10"-12" (25CM - 30CM) | 7" (18CM) |

501 XS, S AVAILABLE IN BLACK ONLY



- Padded pressure sensitive straps are adjustable and comfortable
- Deep pocket design
- Two adjustment points for each strap allows for better fit
- Comfortable fit
- Open end with thumb loop adds comfort and eliminates slippage
- 501 PED constructed of cotton polyester in a fun, colorful, kid friendly print

ORDER No.

DESCRIPTION

661

CLAVICLE SPLINT

LCODE L3650 SUGGESTED


SIZING:

MEASURE FROM CENTER OF SHOULDER BLADES AROUND SHOULDER

MIN

MAX

|    |            |            |
|----|------------|------------|
| XS | 15" (38CM) | 21" (53CM) |
| S  | 19" (48CM) | 28" (71CM) |



- Ideal for clavicle fracture splinting and posture control
- Comfortable foam construction provides proper immobilization
- Contoured urethane foam is hook & loop sensitive allowing maximum adjustability
- Two buckles on rear felt pad are used for tension adjustment

ORDER No.

DESCRIPTION

6674

SHIELDS II® BRACE

LCODE L1820 SUGGESTED

6675

HINGED SHIELDS II® BRACE

LCODE L1820 SUGGESTED

SIZING:

MEASURE 4" ABOVE MID-PATELLA

XS

10" - 11.5" (25CM - 29CM)

S

12" - 13.5" (30CM - 34CM)

MEASURE AT MID-PATELLA

XS

10" - 11" (25CM - 28CM)

S

11.5" - 13" (29CM - 33CM)



- Provides dynamic control of the patella reducing discomfort and enhancing recovery
- Trakaderm® U-shaped buttress allows for treatment of multiple patella indications and fine tuning the position of the patella no matter the malalignment
- One brace fits both left and right for less inventory
- Independent calf and thigh straps adjust for contours of the leg
- Low Profile and extremely durable

ORDER No.

DESCRIPTION

211

12" THREE PANEL KNEE IMMOBILIZER

LCODE L1830 SUGGESTED

SIZING:

FITS MID-PATELLA CIRC.

12" - 24" (30CM - 61CM)



- Provides rigid support for the knee
- 3/8" (9.5mm) pile foam laminate w/ black flannel interior ensures patient comfort
- Two detachable panels for circumferential adjustments for custom fitting
- Posterior and medial/lateral stays are malleable
- Optional patella strap provides extra knee control

ORDER No.

DESCRIPTION

314-BLK

WEBLY ZAP® ANKLE ORTHOSIS

LCODE L1902 SUGGESTED

308A

OPTIONAL STAYS (PAIR)

SIZING:

FITTING BY SHOE SIZE

MEN

NARROW

REGULAR

WIDE

|     |     |         |   |
|-----|-----|---------|---|
| XXS | 6.5 | 5.5 - 6 | 5 |
| XS  | 7.5 | 6.5 - 7 | 6 |
| S   | 8.5 | 7.5 - 8 | 7 |

WOMEN

NARROW

REGULAR

WIDE

|     |      |           |     |
|-----|------|-----------|-----|
| XXS | 7.5  | 6.5 - 7   | 6   |
| XS  | 9    | 7.5 - 8.5 | 7   |
| S   | 10.5 | 8 - 10    | 7.5 |



- Limits ankle inversion/eversion while not restricting plantar and dorsiflexion
- Figure 8 stabilizing straps replicate ankle taping
- Low profile brace fits in shoes
- Useful prophylactically and for ankle sprains
- Seamless design maximizes comfort
- Elastic strap closure improves compression and control of syndesmosis
- Universal fits right or left
- For additional support add stays

ORDER No.

DESCRIPTION

318

RAPID ZAP™ ANKLE ORTHOSIS

LCODE L1902 SUGGESTED

308A

OPTIONAL STAYS (PAIR)

SIZING:

FITTING BY SHOE SIZE

MEN

NARROW

REGULAR

WIDE

|     |     |         |   |
|-----|-----|---------|---|
| XXS | 6.5 | 5.5 - 6 | 5 |
| XS  | 7.5 | 6.5 - 7 | 6 |
| S   | 8.5 | 7.5 - 8 | 7 |

WOMEN

NARROW

REGULAR

WIDE

|     |      |           |     |
|-----|------|-----------|-----|
| XXS | 7.5  | 6.5 - 7   | 6   |
| XS  | 9    | 7.5 - 8.5 | 7   |
| S   | 10.5 | 8 - 10    | 7.5 |



- Unique lacing design allows independent fore foot, mid foot and ankle adjustments
- Figure 8 Straps mimic ankle taping job in a preferred neutral position providing resistance to "roll-over" ankle injury while not restricting walking and running
- Sturdy construction and seamless design offers comfort, application, quality and durability
- Elastic cuff closure secures laces and straps while offering additional compression
- Adjustable lacing panel enables the opening to accommodate wide or swollen ankles
- A brace from acute to sports

39

40

## CUSTOM LOGO

### GIVE YOUR MEDICAL PRACTICE A VISUAL IDENTITY

With the Hely & Weber Custom Logo Program, your practice name will be prominently displayed on any brace you order. Logos may be designed with up to 5 different colors, including your web address and phone number.

We accept camera ready artwork or image files (jpg, gif, ai, tif). Pantone color numbers will expedite order processing.

Please contact your sales representative or customer service for program details and pricing.



## CUSTOM BRACE PROGRAM

### SPECIAL BRACING FOR UNIQUE NEEDS

Our Custom Brace Program is available for all Hely & Weber neoprene knee, elbow, wrist, thumb, ankle and back braces.

Please contact your sales representative or customer service department for program details and pricing. Let us help you meet your patient's specific individual needs.

## GUARANTEE

- Hely & Weber products are 100% guaranteed, and will be replaced in the event of defect in materials or workmanship
- At the discretion of Hely & Weber, products deemed to have been misused, excessively worn or discontinued may not be credited or replaced
- Hely & Weber offers no other guarantee, expressed or implied, including but not limited to the implied guarantees or merchantability and fitness of purpose

## COMMERCIAL CUSTOMERS

### SALES & ORDERING TERMS AND POLICIES

#### ORDER PLACEMENT:

Orders may be called to: 800-654-3241 or 800-221-5465  
Orders may be faxed to: 800-559-5975 or 805-933-0535  
Orders may be mailed to: Hely & Weber, P.O. Box 832,  
Santa Paula, CA 93061-0832

Orders may be e-mailed to: [order@hely-weber.com](mailto:order@hely-weber.com)

#### CUSTOM OR SPECIAL SIZES:

There is a \$40.00 custom fee for non-stocked items  
There is a \$10.00 size fee for XXXL

#### SHIPPING:

Orders are shipped pre-paid via FedEx, UPS, and Priority Mail (USPS) with freight and any handling charges added to the invoice

All orders are shipped F.O.B.  
Santa Paula, CA 93060 and Dallas Fort-Worth, TX 75261

### RETURNS/ EXCHANGES

- All returns/exchanges must be pre-authorized and have a Return Materials Authorization Number printed on the exterior of the box
- Please include any paperwork sent by Hely & Weber with regard to the return/exchange
- All returns are subject to a restocking fee at the discretion of Hely & Weber
- Special orders or custom products may not be returned
- Logoed products may not be returned

### TERMS

- Net 30 days
- 2% charge per month on past due accounts
- We reserve the right to hold shipments on delinquent accounts
- We accept Visa / MasterCard / American Express (no discount given)
- Hely & Weber reserves the right to discontinue products or to change the specifications, designs, prices or the terms of conditions of the sale at any time
- Payments should be mailed to:  
Hely & Weber, P.O. Box 832  
Santa Paula, CA 93061-0832