
Stays for Ankle Braces


Warnings and Instructions: Review carefully, proper application is required

⚠**Warning:** This device will not prevent or eliminate risk of injury. Do Not Overtighten. If swelling, pain, skin irritation, or an unusual reaction occurs, discontinue use immediately and consult your medical professional.

Care: Hand wash using mild soap. Rinse thoroughly. Air dry only. Do not tumble dry. For single patient use only.



1. Secure laces. Center the cutout over the ankle bone and apply hook side of stays to brace.



2. Apply straps over the stays according to directions in instructions for ankle brace.

HELY&WEBER
2018 Hely & Weber


Stays for Ankle Braces

308


Warnings and Instructions: Review carefully, proper application is required

⚠**Warning:** This device will not prevent or eliminate risk of injury. Do Not Overtighten. If swelling, pain, skin irritation, or an unusual reaction occurs, discontinue use immediately and consult your medical professional.

Care: Hand wash using mild soap. Rinse thoroughly. Air dry only. Do not tumble dry. For single patient use only.



1. Secure laces. Center the cutout over the ankle bone and apply hook side of stays to brace.



2. Apply straps over the stays according to directions in instructions for ankle brace.

HELY&WEBER
2018 Hely & Weber

Rev D 10/26